

Deep-Well and 96-Well Plates

Deep-well plates


BRAND deep well plates are excellent for sample storage, high throughput screening (HTS) assays requiring a motherplate, cell and tissue culture, and more. Choose polypropylene for chemical resistance to phenols, chloroform and DMSO, or polystyrene for clarity or surface properties. Plates are compatible with most robotic samplers and automated liquid handling systems

- Polypropylene plates are autoclavable at 250°F (121°C) for 20 minutes, and withstand temperatures as low as -112°F (-80°C).
- Round bottom wells for optimal sample mixing and recovery.
- Wells are alphanumerically coded, with a corner cutaway for easy orientation and well identification.
- Compatible with multichannel pipettes and automated liquid handling systems from all leading manufacturers.
- Plates come in 96-well and 384-well varieties, and are stackable for easy storage.

96-well plates

Polystyrene (PS) 96-well plates are an economical choice for cell and tissue culture, sample storage, immunological assays, and other applications.

- Clear polystyrene plates are available in your choice of flat, round, or conical bottoms.
- Plates have raised rims to simplify sealing with adhesive films, and recessed interwell spaces to reduce well cross-contamination. Optional lids have condensation rings.




BRAND 96-well Polystyrene (L) and Polypropylene (R) Microtitration Plate

ORDERING INFORMATION

Description	Well Volume	Qty./Pack	Cat. No.	2009 List Price
Deep Well Plates				
96-well, polypropylene	0.5mL	48	701346	\$315.00
96-well, polypropylene	1.1mL	24	701350	120.00
96-well, polystyrene	1.1mL	32	701352	152.00
96-well, polypropylene	2.2mL	24	701354	173.00
384-well, polypropylene	0.3mL	48	701355	515.00
Mat Covers				
Mat for 1.1mL volume deep well plates		24	701360	97.00
Mat for 2.2mL volume deep well plates		24	701362	76.00
Sealing Films				
Aluminum, roll		100	781380	241.00
Aluminum, RNase-, DNase-free		100 sheets	781381	194.00
Polypropylene, cross-cut, DMSO-resistant		100 sheets	701366	280.00
Rayon, gas-permeable		100 sheets	701364	191.00
Rayon, gas-permeable, sterile		50 sheets	701365	127.00

BRAND 96-well round well (1.1mL) and square-well (2.2mL) Deep Well Plates, PP



ORDERING INFORMATION

Description	Qty./Pack	Packaging	Cat. No.	2009 List Price
Clear, Flat-bottom PS Plates, 400µL well volume				
Non-sterile	100	20 bags of 5 plates	701304	\$116.00
Clear, Round-bottom PS Plates, 330µL well volume				
Non-sterile	100	20 bags of 5 plates	701300	116.00
Clear, Conical (V)-bottom PS Plates, 330µL well volume				
Non-sterile	100	20 bags of 5 plates	701302	116.00
Lids for Plates Listed Above				
Clear polystyrene plate lids, non-sterile	100	20 bags of 5 lids	701306	104.00