

BRAND Microcentrifuge Tubes

BRAND microcentrifuge tubes fit most common rotors and may be spun up to 20,000 RCF (10,000 RCF for 0.5mL). Manufactured from virgin polypropylene, these highly translucent tubes are free from flash, air bubbles, and streaks. Caps are flat with consistent thickness for trouble-free piercing.

- Tubes have uniform wall thickness for consistent, repeatable results.
- Caps form a tight, leak-free seal that is easy to open, even when wearing gloves.
- BIO-CERT® endotoxin-, RNase-, ATP-, and DNA-free microcentrifuge tubes are available for critical applications.
- Colored tubes simplify sample identification.
- Centrifuge Tube Mini Coolers hold twelve microtubes from 0.5mL to 2.0mL in size and, when pre-frozen, maintain the rated temperature for approximately an hour (45 minutes for -70°C model).



BRAND 1.5mL Microcentrifuge tubes in colors



0.5mL Microcentrifuge tubes



1.5mL Microcentrifuge tubes



2.0mL Microcentrifuge tubes



128-tube Rack



Centrifuge Tube Mini Cooler, 0°C



20-tube Rack

ORDERING INFORMATION

Tube Vol mL	Color	OD x H mm	Qty./Pack	Packaging	Cat. No.	2011 List Price
Non-sterile Microcentrifuge Tubes with Lids						
0.5	Clear	8 x 32	1000	1 pack of 1000 tubes	780507	\$61.40
1.5	Clear	11 x 41	500	1 pack of 500 tubes	780500	12.40
1.5	Yellow	11 x 41	500	1 pack of 500 tubes	780521	12.40
1.5	Blue	11 x 41	500	1 pack of 500 tubes	780522	12.40
1.5	Green	11 x 41	500	1 pack of 500 tubes	780523	12.40
1.5	Orange	11 x 41	500	1 pack of 500 tubes	780524	12.40
1.5	Amber	11 x 41	500	1 pack of 500 tubes	780525	12.40
2.0	Clear	11 x 40	500	1 pack of 500 tubes	780550	47.40

BIO-CERT® RNase-, DNase, DNA-free Microcentrifuge Tubes with Lids

1.5	Clear	11 x 40	450	30 15-tube blister packs	780400	113.00
-----	-------	---------	-----	--------------------------	--------	--------

Non-sterile Microcentrifuge Tubes without Lid

0.5	Clear	11 x 39.5	1000	1 pack of 1000 tubes	780505	48.00
-----	-------	-----------	------	----------------------	--------	-------

Accessories

128-tube Rack	5	1 pack of 5 racks	4341052	90.00
84-tube Rack	5	1 pack of 5 racks	4341002	90.00
20-tube Rack	1	Each	780605	16.40
0.5mL Adapters	20	1 pack of 20	780608	21.60
Centrifuge Tube Mini Cooler, 0°C	1	Each	114930	111.00
Centrifuge Tube Mini Cooler -20°C	1	Each	114935	111.00
Centrifuge Tube Mini Cooler, -70°C	1	Each	114940	111.00