

# BRAND Microcentrifuge Tubes

New BRAND microcentrifuge tubes fit most common rotors and may be spun up to 20,000 RCF (10,000 RCF for 0.5mL). Manufactured from virgin polypropylene, these highly translucent tubes are free from flash, air bubbles, and streaks. Caps are flat with consistent thickness for trouble-free piercing.

- Tubes have uniform wall thickness for consistent, repeatable results.
- Caps form a tight, leak-free seal that is easy to open, even when wearing gloves.
- BIO-CERT® endotoxin-, RNase-, ATP-, and DNA-free microcentrifuge tubes are available for critical applications.
- Colored tubes simplify sample identification.
- Centrifuge Tube Mini Coolers hold twelve microtubes from 0.5mL to 2.0mL in size and, when pre-frozen, maintain the rated temperature for approximately an hour (45 minutes for -70°C model).



BRAND microcentrifuge tubes in colors

## ORDERING INFORMATION

Tube Vol mL	Color	OD x H mm	Qty./ Pack	Packaging	Cat. No.	2009 List Price
<b>Non-sterile Microcentrifuge Tubes with Lids</b>						
0.5	Clear	8 x 32	1,000	1 pack of 1,000 tubes	780507	\$60.00
1.5	Clear	11 x 41	500	1 pack of 500 tubes	780500	12.00
1.5	Yellow	11 x 41	500	1 pack of 500 tubes	780521	12.00
1.5	Blue	11 x 41	500	1 pack of 500 tubes	780522	12.00
1.5	Green	11 x 41	500	1 pack of 500 tubes	780523	12.00
1.5	Orange	11 x 41	500	1 pack of 500 tubes	780524	12.00
1.5	Amber	11 x 41	500	1 pack of 500 tubes	780525	12.00
2.0	Clear	11 x 40	500	1 pack of 500 tubes	780550	46.00

## BIO-CERT® RNase-, DNase, DNA-free Microcentrifuge Tubes with Lids

1.5	Clear	11 x 40	450	30 15-tube blister packs	780400	35.00
-----	-------	---------	-----	--------------------------	--------	-------

## Non-sterile Microcentrifuge Tubes without Lid

0.5	Clear	11 x 39.5	1,000	1 pack of 1,000 tubes	780505	47.00
-----	-------	-----------	-------	-----------------------	--------	-------

## Accessories

128-tube Rack	5	1 pack of 5 racks	4341052	88.00
84-tube Rack	5	1 pack of 5 racks	4341002	88.00
20-tube Rack	1	Each	780605	16.00
0.5mL Adapters	20	1 pack of 20	780608	21.00
Centrifuge Tube Mini Cooler, 0°C	-	-	-	-
1	Each	114930	106.00	
Centrifuge Tube Mini Cooler -20°C	-	-	-	-
1	Each	114935	106.00	
Centrifuge Tube Mini Cooler, -70°C	-	-	-	-
1	Each	114940	106.00	



1.5mL Microcentrifuge tubes



1.1mL Microcentrifuge tubes



2.0mL Microcentrifuge tubes



128-tube rack



Centrifuge Tube Mini Cooler, 0°C



BRAND