

BRAND PCR* Tubes and Plates


BRAND PCR* tubes and plates are manufactured from high-quality virgin resins for excellent, reproducible results in most PCR* applications. This collection of individual tubes, tube strips, cap strips and multi-well plates are manufactured without lubricants or chemical additives that can contaminate reagents and samples.

- PCR* tubes and plates from BRAND are compatible with most thermal cyclers. Tube strips and plates are available in formats compatible with 24-, 32-, 96- and 384-well thermal blocks.
- BRAND PCR* tubes and plates are autoclavable and RNase-, DNase-, and DNA-free.
- The thin, uniform walls of BRAND PCR* tubes and plates ensure consistent, precise thermal transfer for excellent results.
- All tubes with flat caps, and plates sealed with polyester film (Cat.No. 781391) are suitable for real-time PCR*.
- Individual tubes have attached caps that form a tight seal to help prevent sample evaporation, yet are easy to open and close.

- Innovative TPE mats for 96-well plates are autoclavable, flexible and, when used with BRAND PCR* plates, dramatically reduce sample evaporation. Mat material can be easily pierced with a pipette tip.
- All sealing films are DNA-, RNase- and DNase-free.
- BRAND's PCR* rack securely holds 96-well PCR plates, tube strips, and individual tubes for easy handling, and is alphanumerically marked for sample identification. The PCR* box holds the rack for storage, and withstands temperatures from -80°C to +121°C.

NEW! PCR* Mini Cooler helps protect samples from warming. Insulating gel maintains 4°C temperature for approximately 3 hours, and changes color from violet to pink when temperature reaches 7°C. Suitable for 0.2mL tubes, strips and plates. Transparent lid included.

* The Polymerase Chain Reaction (PCR) is a patented process of Hoffman-La Roche. Use of the PCR process requires a license.



Uniform Thickness Provides Reliable, Reproducible Results.





96-well PCR* plate, unskirted with TPE mat.



PCR* Rack and Box



BRAND PCR* tubes



BRAND PCR* tube strips



PCR* Mini-Cooler

ORDERING INFORMATION

Description	Qty./Pack	Cat. No.	2009 List Price
PCR* Tubes			
0.2mL tubes with domed cap, clear	1000	781300	\$93.00
0.2mL tubes with flat cap, clear	1000	781305	49.00
0.2mL tubes with flat cap, rose	1000	781301	49.00
0.2mL tubes with flat cap, yellow	1000	781302	49.00
0.2mL tubes with flat cap, green	1000	781303	49.00
0.2mL tubes with flat cap, blue	1000	781304	49.00
0.5mL tubes with flat cap, clear	1000	781310	49.00
0.5mL tubes with flat cap, rose	1000	781311	49.00
0.5mL tubes with flat cap, yellow	1000	781312	49.00
0.5mL tubes with flat cap, green	1000	781313	49.00
0.5mL tubes with flat cap, blue	1000	781314	49.00
PCR* 8-Strip Tubes			
NEW! 8 x 0.2mL tube strip with individual attached flat caps			
	120	781332	190.00
NEW! 8 x 0.15mL low profile tube strip with individual attached flat caps			
	120	781333	218.00
NEW! 8 x 0.2mL tube strip with attached domed strip cap			
	125	781330	189.00
8 x 0.2mL tube strip, clear	125	781320	63.00
8 x 0.2mL tube strip, rose	125	781321	63.00
8 x 0.2mL tube strip, yellow	125	781322	63.00
8 x 0.2mL tube strip, green	125	781323	63.00
8 x 0.2mL tube strip, blue	125	781324	63.00
PCR* 8-Strip Caps			
NEW! 8-strip flat caps, clear	125	781334	51.00
8-strip domed caps, clear	125	781340	31.00
8-strip domed caps, rose	125	781341	31.00
8-strip domed caps, yellow	125	781342	31.00
8-strip domed caps, green	125	781343	31.00
8-strip domed caps, blue	125	781344	31.00

ORDERING INFORMATION

Description	Qty./Pack	Cat. No.	2009 List Price
PCR* 12-Strip Tubes			
12 x 0.2mL tube strip, clear	125	781280	\$180.00
12 x 0.2mL tube strip, rose	125	781281	180.00
12 x 0.2mL tube strip, yellow	125	781282	180.00
12 x 0.2mL tube strip, green	125	781283	180.00
12 x 0.2mL tube strip, blue	125	781284	180.00
PCR* 12-Strip Caps			
12-strip domed caps, clear	125 strips	781290	45.00
12-strip domed caps, rose	125 strips	781291	45.00
12-strip domed caps, yellow	125 strips	781292	45.00
12-strip domed caps, green	125 strips	781293	45.00
12-strip domed caps, blue	125 strips	781294	45.00
PCR* Plates			
96-well plate, non-skirted	5 bags of 10	781350	163.00
96-well plate, semi-skirted	5 bags of 10	781400	213.00
384-well plate, skirted	5 bags of 10	781345	260.00
PCR* Mats and Sealing Films			
TPE-mat for BRAND 96-well PCR* plates	5	781405	46.00
NEW! Highly-transparent polyester film for real-time PCR*			
	100 sheets	781391	260.00
Polypropylene film	100 sheets	781390	385.00
Aluminum film	Roll of 100 sheets	781380	241.00
NEW! Aluminum film, RNase-, DNase-free	100 sheets	781381	194.00
Racking and Storage			
PCR*-Rack	5	781370	38.00
PCR*-Box	5	781360	63.00
NEW! PCR* Mini-Cooler	2	781260	113.00

* The Polymerase Chain Reaction ("PCR") is covered by patents owned by Hoffmann-La Roche (see trademark index).