

PD-Tip™ Syringe Tips

BRAND encoded PD-Tip™ syringe tips are manufactured from high-quality virgin plastics for use with all standard repeating pipettes in a wide variety of applications. PD-Tip™ syringe tips are made without lubricants or chemical additives to prevent reagent contamination.

- **Compatible with Most Pipettes:** PD-Tip™ syringe tips are an excellent choice for the BRAND HandyStep® and HandyStep® electronic repeating pipettes; the black-handled Eppendorf® Repeater® and EDOS® 5221*; the Rainin® AutoRep™ E and AutoRep™ M; and all other repeating pipettes using standard tips. PD-Tip™ syringe tips are encoded for automatic recognition by the BRAND HandyStep® electronic repeating pipette and the Rainin® AutoRep™ E pipette.
- **Reliable:** A precision seal between the piston and cylinder ensure smooth operation, increased accuracy, and increased precision. Extended piston design minimizes dead volume and maximizes fluid utilization.
- **Manufactured from High Quality Materials:** PD-Tip™ syringe tips are made entirely from high-quality virgin polyethylene and polypropylene without lubricants or other chemical additives. (Piston of 0.1mL tip is liquid crystal polymer-LCP).

The Perfect
Tips for Your
Repeating
Pipette.



*The Eppendorf Repeater Plus, Repeater Pro, Stream, and XStream are designed to preclude their operation with any tips other than Eppendorf Combitips Plus.

ISO/GLP Calibration Support

BRAND PD-Tip™ syringe tips are supplied with a batch Certificate of Performance which documents tip accuracy and precision testing.

Preserve Sample Integrity with BIO-CERT® Quality Sterile Tips

BIO-CERT® PD-Tip™ syringe tips are sterile, individually wrapped and certified free of DNA, RNase, endotoxins, and ATP to meet the most stringent sample integrity requirements. For more information on the sterility and purity levels, see page 19.

ORDERING INFORMATION

Capacity, mL	Cat. No.	2009 List Price
PD-Tip™ Syringe Tips, non-sterile		
0.1mL, pack of 100	702402	\$110.00
0.5mL, pack of 100	702370	103.00
1.25mL, pack of 100	702372	103.00
2.5mL, pack of 100	702374	103.00
5mL, pack of 100	702376	103.00
12.5mL, pack of 100	702378	103.00
25mL, pack of 50 (one adapter included)	702380	64.00
50mL, pack of 25 (one adapter included)	702382	36.00
PD-Tip™ Syringe Tips, non-sterile, assorted pack, 20 each of 0.5, 1.25, 2.5, 5 and 12.5mL sizes	702368	103.00
Adapter for 25mL and 50mL PD-Tip Syringe Tips, non-sterile, PP, autoclavable, pack of 10	702398	83.00
PD-Tip™ Syringe Tips, BIO-CERT® sterile, individually wrapped		
0.1mL, pack of 100	702683	147.00
0.5mL, pack of 100	702684	140.00
1.25mL, pack of 100	702686	140.00
2.5mL, pack of 100	702688	140.00
5mL, pack of 100	702690	140.00
12.5mL, pack of 100	702692	140.00
25mL, pack of 25 (one adapter included)	702694	40.00
50mL, pack of 25 (one adapter included)	702696	41.00
Adapter for 25 and 50mL PD-Tip™ Syringe Tips, BIO-CERT® sterile, PP, autoclavable, pack of 5	702699	64.00