

# QuikSip™

QuikSip™ bottletop aspirator

The QuikSip™ is a convenient bottletop aspirator that removes supernatants, cell culture media, and other biological liquids from dishes, flasks, culture bottles, centrifugation tubes, multiwell plates, and more. It is small, portable, and doesn't require an external vacuum or power source.

- **Moves Easily Between Labs:** The QuikSip™ works without any external vacuum source, allowing it to be easily transferred between workstations. Each plunger stroke powers up to 25mL of aspiration.
- **Provides Fingertip Aspiration Control:** The QuikSip™ includes a unique thumb-operated controller that is compatible with disposable pipette tips, micropipettes, and glass Pasteur pipettes.
- **Adapts to Sterile Environments:** All wetted parts are autoclavable (except piston/cylinder cartridge). Sterile piston/cylinders are available for use in sterile applications.
- **Fits Most Reagent Containers:** Mounts on most common reagent bottles with its 45mm thread and 45/33mm and 45/38mm adapters (included). Other thread adapters available.

## The Economy Aspirator for Biological Fluids.



## Technical Data

### Operating limits

- Liquids unsuitable for use with the QuikSip™ bottletop aspirator include, but are not limited to:
  - Concentrated or oxidizing acids;
  - Non-polar solvents, such as hydrocarbons and halogenated hydrocarbons;
  - Hydrofluoric acid (HF);
  - Those that attack PP, PE, EPDM and silicone.
- Temperature of instrument and reagents must be between 59° and 104°F (15° and 40°C).
- Piston/cylinder cartridge and pump housing are not autoclavable.

### Items supplied

QuikSip™ aspirator includes:

- QuikSip™ controller with 78in. (200cm) of silicone tubing
- Spare piston/cylinder cartridge
- PP thread adapters (GL 45/33 and GL 45/38)
- Operating manual
- One-year warranty



QuikSip™ aspirator shown with optional autoclavable 8-channel manifold.

### ORDERING INFORMATION

Description	Cat. No.	2009
		List Price
QuikSip™ Bottletop Aspirator	4724150	\$285.00
8-channel manifold; PP, autoclavable, pack of 1	704526	13.00
25mL replacement piston/cylinder, cartridge, pack of 3	704504	56.00