



**Improved!**

# Stainless Steel Support

Stainless steel support jacks from BrandTech Scientific are rugged supports that are excellent for a wide variety of critical laboratory applications. Constructed entirely from stainless steel with POM rollers and PF31 phenolic knob, these support jacks are exceptionally stable and well-suited for corrosive environments such as fume hoods, water baths, marine and industrial labs, and more.

- **Supports a Wide Variety of Items:** Stainless steel support jacks are available in a range of sizes and load capacities up to 176 pounds for most lab applications.
- **Improved Lift Mechanism:** Provides increased stability and torsional rigidity.
- **Helps Reduce Contamination:** Support jacks feature smooth surfaces that are easy to clean and sanitize. All support jacks are fully autoclavable at 121°C for sterile applications.
- **Reaches New Heights:** These stainless steel support jacks are now also available in Maxi-Lift models. Maxi-Lift stainless steel support jacks have a maximum height that is approximately 30% higher than traditional support jacks for special applications.
- **Offers a Wide Variety of Options:** BrandTech Scientific offers support jacks manufactured with anodized aluminum and powder-coated aluminum decks (on stainless steel lift mechanisms) for less critical applications.
- **Two-year warranty on all support jacks**



## Rugged Laboratory Supports for Critical Environments.

# Jacks

## Optional Handle

An optional ratcheting handle is available to simplify use of BrandTech's larger support jacks.

This handle simply slides onto the handwheel of 30 x 30cm and 40 x 40cm support jacks to greatly reduce the effort required for raising and lowering the instrument.



## Also Available

Support jacks with epoxy powder-coated aluminum and anodized aluminum decks, with stainless steel lift mechanisms.

## Epoxy Powder Coat



## Anodized Aluminum



## ORDERING INFORMATION

Description Plate, L x W, cm	Minimum Height, cm	Maximum Height, cm	Maximum Static Load, kg	Cat. No.	2011 List Price
<b>Stainless Steel Support Jacks</b>					
10 x 10	5.5	12	10	B11115	\$218.00
16 x 13	6.0	27.5	20	B11120	233.00
20 x 20	6.0	27.5	20	B11130	288.00
24 x 24	6.0	27.5	20	B11140	307.00
30 x 30	13.0	47	80	B11180	1,495.00
40 x 40	13.0	47	80	B11190	1,678.00

## Maxi-Lift Stainless Steel Support Jacks

16 x 13	7.5	40	20	B11122	288.00
20 x 20	7.5	40	20	B11132	352.00

## Anodized Aluminum Deck Support Jacks

10 x 10	5.5	12	10	B11016	168.00
16 x 13	6.0	27.5	20	B11021	188.00
20 x 20	6.0	27.5	20	B11031	228.00
24 x 24	6.0	27.5	20	B11041	253.00
30 x 30	13.0	47	80	B11081	1,086.00
40 x 40	13.0	47	80	B11091	1,282.00

## Powder-coated Aluminum Deck Support Jacks

10 x 10	5.5	12	10	B11015	154.00
16 x 13	6.0	27.5	20	B11020	174.00
20 x 20	6.0	27.5	20	B11030	209.00
24 x 24	6.0	27.5	20	B11040	231.00
30 x 30	13.0	47	80	B11080	952.00
40 x 40	13.0	47	80	B11090	1,174.00

## Accessories

### Optional Ratchet Wrench

for B11180, B11190, B11081, B11091, B11080, B11090	B11089	57.00
--	--------	-------

1cm = 0.394 inches, 1kg = 2.21 pounds