

NEW! Test Tube Racks

Stainless Steel Test Tube Racks

Racks are constructed entirely of rugged, corrosion-resistant 304 stainless steel. Highly polished surfaces reduce porosity and reinforce the corrosion-resistant properties.

ORDERING INFORMATION

Max.							
Tube Diam. in. (cm)	Max. # of Tubes	Format	Length in. (cm)	Width in. (cm)	Height in. (cm)	Cat. No.	2009 List Price
Test Tube Racks (each)							
0.7 (1.7)	12	2x6	5.7 (14.5)	2.0 (5.2)	3.5 (9.0)	B8210	\$49.00
0.7 (1.7)	24	2x12	10.9 (27.8)	2.0 (5.2)	3.5 (9.0)	B8212	53.00
0.7 (1.7)	48	4x12	11.3 (28.7)	3.8 (9.7)	3.5 (9.0)	B8214	93.00
0.8 (2.0)	12	2x6	6.5 (16.4)	2.4 (6.0)	3.5 (9.0)	B8220	63.00
0.8 (2.0)	24	2x12	12.2 (31.0)	2.4 (6.0)	3.5 (9.0)	B8222	81.00



Test Tube Racks

ORDERING INFORMATION

Max.							
Tube Diam. in. (cm)	Max. # of Tubes	Format	Length in. (cm)	Width in. (cm)	Height in. (cm)	Cat. No.	2009 List Price
Test Tube Racks, Z-style (each)							
0.5 (1.3)	50	10x5	2.0 (5.0)	6.5 (16.5)	0.5 (1.3)	B8290	\$42.00
0.7 (1.7)	50	10x5	2.0 (5.0)	8.7 (22.0)	0.7 (1.7)	B8291	42.00
0.5 (1.3)	25	5x5	2.0 (5.0)	3.3 (8.5)	0.5 (1.3)	B8292	31.00
0.7 (1.7)	25	5x5	2.0 (5.0)	4.3 (11.0)	0.7 (1.7)	B8293	31.00



Test Tube Racks, Z-style

Polypropylene Test Tube Racks

Sturdy polypropylene racks with alphanumeric positions for easy sample identification. Racks are stackable, non-conductive and will not float in waterbaths. They are suitable for use in temperatures from -20°C to 90°C. Racks are supplied in packs of five. They are supplied in three pieces for easy and permanent assembly.

ORDERING INFORMATION

Max.										
Tube Diam. in. (mm)	Format	Length in. (mm)	Width in. (mm)	Height in. (mm)	White Cat. No.	Blue Cat. No.	Red Cat. No.	Yellow Cat. No.	2009 List Price	
Test Tube Racks, PP (pack of 5)										
13	6x14	265	126	75	4340000	4340001	4340002	4340003	\$68.00	
16	5x11	265	126	75	4340060	4340061	4340062	4340063	68.00	
18	5x11	265	126	75	4340010	4340011	4340012	4340013	68.00	
20	4x10	265	126	75	4340020	4340021	4340022	4340023	68.00	
25	4x8	265	126	88	4340030	4340031	4340032	4340033	68.00	
30	3x7	265	126	88	4340040	4340041	4340042	4340043	68.00	



Test Tube Racks, PP

