


NEW! Pipette Tips

BRAND pipette tips are made from high purity virgin polypropylene and are designed for use in a wide variety of pipetting applications. They are tested for compatibility with BRAND, Eppendorf®, Gilson®, Rainin®. (See individual tip descriptions on page 15 for specifics). Suitable for use with instruments from many other manufacturers, except LTS®.

- **Oleamide- and DiHEMA-FREE!** BRAND tips are now free of these commonly used plastic additives which can affect biological testing, and cause false results.
- **Right for Most Applications:** BRAND pipette tips are available in a variety of sizes and configurations for delivering volumes from 0.1 µL to 10mL. Tips have clean, hydrophobic surfaces to reduce sample retention and increase accuracy and reproducibility. Tip orifice is precisely centered for directional accuracy.
- **Made from High-Quality Materials:** Tips are precision molded from virgin polypropylene, without lubricants that can contaminate samples. Pigments, when used, are cadmium-free to preserve sample integrity.
- **Supplied in Easy-To-Use Packaging:** BRAND tips are available in a variety of packaging options to accommodate the requirements of any application.



**Superior Quality Tips.
Affordably Priced for Value.**



Quality Tips

Type of Tip

[A] 0.1-20 μ L, nano-cap™ tips (neutral)

Compatible with the 0.1-1 μ L Transferpette® pipette as well as most pipettes using ultra-micro tips. Ideal for molecular biological applications, it features a capillary rise that is visible to the naked eye at 0.1 μ L!



[B] 0.5-20 μ L, ultra-micro tips (neutral)

The elongated, slender design of these tips reduces accidental contact with plate and tube surfaces. Graduations at 2 μ L and 10 μ L.



[C] NEW! 2-200 μ L tips (yellow or neutral)

Newly redesigned, these tips now use 25% less plastic than their predecessors and have graduations, improved fit and translucency. Fits most major brands of pipettes. Graduations at 10 μ L and 100 μ L.



[D] 5-300 μ L tips (neutral)

Designed for plate washing and reverse pipetting. Fit most 200 μ L pipettes, including multichannels.



[E] NEW! 50-1000 μ L tips (blue or neutral)

Newly redesigned tips are optimized for best fit with leading pipettes, and use 25% less plastic. Graduations at 250 μ L, 500 μ L and 1000 μ L give a visual indication of volume.



[F] 50-1250 μ L tips (neutral)

Extended tip length (82mm) and slender shape allows pipetting into narrow vessels. Molded graduations at 250, 500 and 1000 μ L permit easy verification of approximate volumes. Also suitable for reverse pipetting at higher volumes. Fit most instruments using 1000 μ L tips.



[G] 0.5-5mL tips (neutral)

Slim design (160mm long; only 9mm diameter) allows for pipetting into narrow vessels. Tested for use with BRAND and Finnpipette® pipettes.



[H] 1-10mL tips (neutral)

Ideal for working with the Transferpette® S 10mL. Compatible with Eppendorf® and Gilson® systems.




BRAND pipette tips

Innovative Packaging Options

BRAND pipette tips are packaged bulk or racked. Racked tips are available in boxes, refill insert packs, and as insert stacks.

All rack-packed pipette tips are made from neutral polypropylene. For easy identification, the tip insert trays are colored to indicate tip size.

Identical dimensions of 60- and 96-tip insert trays allow both the non-sterile and BIO-CERT® Tip-Refills to fit the Tip-Box SL, the Tip-Box PP, and the Tip-Box PC.



Tip-Box SL

Rugged polycarbonate tip box with an innovative lid that slides open left or right, forward or back for convenient eight- or twelve-channel pipetting. Keeps unused tips covered and reduces the risk of contamination. Autoclavable at 121°C (250°F) for 20 minutes. Fits Tip-Refill-60 and Tip-Refill-96.



Tip-Box PP

This space-saving polypropylene tip box features a push-on lid and is designed for long life, even when subjected to repeated autoclavings at 121°C (250°F). Tall version refillable with 60- and 96-tip inserts. Short version refillable with 96-tip Tip-Refill only.



Tip-Box-100, PP

This new box size provides convenient rack storage of 1000µL and 1250µL pipette tips in racks of 100. Refillable only with Tip-Refill-100.

ORDERING INFORMATION

Description	Cat. No.	2011 List Price
Tip-Box SL, empty	702200	\$46.00



Tip-Box, 5mL, 10mL

Designed specifically for BRAND 5mL and 10mL tips. Withstands repeated autoclaving at 121°C (250°F). Holds 28 or 18 tips, respectively.



Tip-Refill, Tip-Refill-100

Refill inserts for BRAND tip boxes in space-saving dispenser cartons. Tip-Refill-100 fits Tip-Box-100 only.



Tip-Stack™

Tip-Stack™ is a space-saving tip-refill system that reduces packaging waste, and simplifies refilling the Tip-Box PC, Tip-Box SL and the Tip-Box PP. It consists of ten stacked and nested 96-tip inserts filled with tips, and one Tip-Box PP. Newly available in ultra-micro, 200µL and 1000µL sizes.

Pipette Tips: Ordering Information

ORDERING INFORMATION

Capacity/ Volume	Packaging	# of tips	Description	Cat. No.	2011 List Price
0.1-20µL nano-cap™ tips (dark gray/neutral)*					
0.1-20µL	Bulk	2000	2 bags of 1000 tips	702504	\$83.00
0.1-20µL	Tip-Box SL	96	1 polycarbonate tip box with sliding lid	702202	49.00
0.1-20µL	Tip-Box PP	480	5 polypropylene tip boxes filled with 96 tips each	702408	72.00
0.1-20µL	Tip-Refill	960	10 filled 96-tip inserts	702312	94.60
0.5-20µL ultra-micro tips (dark gray/neutral)*					
0.5-20µL	Bulk	2000	2 bags of 1000 tips	702526	63.00
0.5-20µL	Bulk	10000	10 bags of 1000 tips	702565	278.00
0.5-20µL	Tip-Box SL	96	1 polycarbonate tip box with sliding lid	702204	49.00
0.5-20µL	Tip-Box PP	480	5 polypropylene tip boxes filled with 96 tips each	702411	50.00
0.5-20µL	Tip-Stack™	960	10 filled and stacked 96-tip inserts, and one polypropylene box	712402	82.00
0.5-20µL	Tip-Refill	960	10 filled 96-tip inserts	702350	82.00
NEW! 2-200µL tips (yellow/neutral)*					
2-200µL	Bulk	1000	1 bag of 1000 tips	712516	16.00
2-200µL	Bulk	10000	10 bags of 1000 tips	712556	125.00
2-200µL	Tip-Box SL	96	1 polycarbonate tip box with sliding lid	712206	47.00
2-200µL	Tip-Box PP	480	5 polypropylene tip boxes filled with 96 tips each	712413	36.00
2-200µL	Tip-Refill	960	10 filled 96-tip inserts	712315	57.60
2-200µL	Tip-Stack™	960	10 filled and stacked 96-tip inserts, and one polypropylene box	712405	58.60
5-300µL tips (green/neutral)*					
5-300µL	Bulk	1000	1 bag of 1000 tips	702533	26.00
5-300µL	Tip-Box SL	96	1 polycarbonate tip box with sliding lid	702208	47.00
5-300µL	Tip-Box PP	480	5 polypropylene tip boxes filled with 96 tips each	702415	47.00
5-300µL	Tip-Refill	960	10 filled 96-tip inserts	702353	75.00
NEW! 50-1000µL tips (blue/neutral)*					
50-1000µL	Bulk	1000	2 bags of 500 tips	712521	18.40
50-1000µL	Bulk	5000	10 bags of 500 tips	712561	72.00
50-1000µL	Tip-Box SL	60	1 polycarbonate tip box with sliding lid	712210	48.00
50-1000µL	Tip-Stack™	960	10 filled and stacked 96-tip inserts, and one polypropylene box	712408	93.60
50-1000µL	Tip-Refill	960	16 filled 60-tip inserts	712320	79.00
50-1000µL	Tip-Box-100	500	5 filled 100 tip boxes	712414	46.00
50-1000µL	Tip-Refill-100	1000	10 filled 100-tip inserts	712322	93.60
50-1250µL tips (orange/neutral)*					
50-1250µL	Bulk	1000	2 bags of 500 tips	702581	85.00
50-1250µL	Tip-Box-100	500	5 filled 100-tip boxes	702418	58.60
50-1250µL	Tip-Refill-100	1000	10 filled 100-tip inserts	702324	120.00
0.5-5mL (500-5000µL) tips (neutral)*					
0.5-5mL	Bulk	200	1 bag of 200 tips	702595	37.00
0.5-5mL	Bulk	1000	5 bags of 200 tips	702600	174.00
0.5-5mL	Tip-Box	28	1 polypropylene box of 28 tips	702605	39.00
1-10mL (1000-10000µL) tips (neutral)*					
1-10mL	Bulk	200	2 bags of 100 tips	702603	74.00
1-10mL	Bulk	1000	10 bags of 100 tips	702604	332.00
1-10mL	Tip-Box	18	1 polypropylene box of 18 tips	702608	77.00

* Bulk 200µL and 1000µL tips feature cadmium-free pigments; other bulk tips are pigment-free. Racked tips have neutral color, pigment-free tips in color-coded carrier plates.

Type of Tip

[A] 0.1-20µL nano-cap™ tips
(dark gray/neutral)*



[B] 0.5-20µL ultra-micro tips
(dark gray/neutral)*



[C] 2-200µL tips
(yellow/neutral)*



[D] 5-300µL tips
(green/neutral)*



[E] 50-1000µL tips
(blue/neutral)*



[F] 50-1250µL tips
(orange/neutral)*



[G] 0.5-5mL (500-5000µL) tips
(neutral)*



[H] 1-10mL (1000-10000µL) tips
(neutral)*




NEW! BIO-CERT® Pipette Tips

BRAND BIO-CERT® pipette tips are sterile pipette tips that are certified free of DNA, RNase, ATP, and pyrogens to preserve sample integrity. They are made from high purity virgin polypropylene and are designed for sensitive applications. They are compatible with both robotic and manual pipetting systems and are tested for compatibility with BRAND, Eppendorf®, Gilson®, Rainin® (except LTS™), and Thermo Labsystems/Finnpipette® pipettes.

- **Suitable for Critical Applications:** All BIO-CERT® pipette tips are certified to be sterile and free of endotoxins, DNA, RNase and ATP.
- **Ideal for Precise, Accurate Pipetting:** Tips are straight and have clean, uniform bores and hydrophobic surfaces to reduce sample retention and increase accuracy and reproducibility.
- **Available in a Variety of Sizes:** BRAND BIO-CERT® tips are available in a variety of sizes to deliver volume from 0.1µL to 1250µL.
- **“Green” Laboratory Packaging:** Refill inserts available in ultra-low-waste, recyclable PET bubble packaging.
- **BIO-CERT® Tip-Box S & Tip-Box 100-S (sterile):** Filled boxes of sterile tips. Tip-Boxes holding 96 racked tips may be refilled with BIO-CERT® Tip-Refills. Tip-Box 100-S holds one hundred 1000µl or 1250µL tips in a disposable box.

BIO-CERT®



**Superior Quality Tips
Certified Pure for Confidence.**



BIO-CERT® Products & Ordering Information

Sensitive applications like PCR, DNA/RNA purification, or DNA sequencing require the highest quality disposable plastics. The advantages of BIO-CERT® pipette tips are:

- **Sterility:** BIO-CERT® products are sterilized according to ISO 11137 and AAMI guidelines with beta-radiation, using a minimum radiation dose of 12.1kGy. This ensures a Sterility Assurance Level (SAL) of 10⁻⁶. The sterility meets USP 29 and Ph. Eur. requirements for sterility.
- **DNA- and RNase-Free:** BIO-CERT® products are DNA-free (less than 4 x 10⁻¹²g/item) to prevent false positive results during PCR applications. These products are free of RNases (less than 8.6 x 10⁻¹⁵g/item) to prevent degradation of RNA during purification processes.
- **Free of Endotoxins:** The endotoxin concentration in BIO-CERT® products is tested by a gel-clotting Limmulus Amoebocyte Lysate test. The detection limit is 0.01 I.E./mL. This gives a documented endotoxin level of less than 1.1 x 10⁻¹²g/item.
- **ATP-Free:** All products are supplied with an ATP concentration of less than 1 x 10⁻¹⁵g/item. BIO-CERT® products are suitable for ATP detection using bioluminescence systems.
- **Lot Tested:** Lots are tested for sterility and pyrogens (endotoxins) and monitored for DNA, RNase, and ATP concentration.
- **Certified:** All products are supplied with a Certificate of Analysis documenting purity test results.
- **Also with Filters:** BRAND also offers BIO-CERT® purity in filter pipette tips. See page 21 for a description and ordering information.



BIO-CERT® Tip-Set (sterile)

Comes complete with one filled sterile Tip-Box S plus four BIO-CERT® tip refill inserts. 1000µL Tip-Set S comes with seven BIO-CERT® tip refill inserts.



BIO-CERT® Tip-Refill (sterile)


Filled tip inserts in low-waste recyclable PET packaging. Designed to refill 60- or 96-tip format BRAND Tip-Boxes.

ORDERING INFORMATION

Capacity/ Volume	Packaging	# of tips	Description: All tips BIO-CERT® pure	Cat. No.	2011 List Price
0.1-20µL nano-cap™ tips (dark gray/neutral)					
0.1-20µL	Tip-Refill S	960	10 filled 96-tip inserts	702428	\$148.00
0.1-20µL	Tip-Set S	480	5 filled 96-tip inserts, and one polypropylene box	702419	94.00
0.1-20µL	Tip-Box S	960	10 polypropylene boxes filled with 96 tips each	702440	148.00
0.5-20µL ultra-micro tips (dark gray/neutral)					
0.5-20µL	Tip-Refill S	960	10 filled 96-tip inserts	702431	130.00
0.5-20µL	Tip-Set S	480	5 filled 96-tip inserts, and one polypropylene box	702421	81.00
0.5-20µL	Tip-Box S	960	10 polypropylene boxes filled with 96 tips each	702442	130.00
NEW! 2-200µL tips (yellow/neutral)					
2-200µL	Tip-Refill S	960	10 filled 96-tip inserts	712433	94.00
2-200µL	Tip-Set S	480	5 filled 96-tip inserts, and one polypropylene box	712423	64.60
2-200µL	Tip-Box S	960	10 polypropylene boxes filled with 96 tips each	712444	99.00
5-300µL tips (green/neutral)					
5-300µL	Tip-Refill S	960	10 filled 96-tip inserts	702435	99.00
5-300µL	Tip-Set S	480	5 filled 96-tip inserts, and one polypropylene box	702425	65.00
5-300µL	Tip-Box S	960	10 polypropylene boxes filled with 96 tips each	702446	99.00
NEW! 50-1000µL tips (blue/neutral)					
50-1000µL	Tip-Box 100-S	1000	10 filled 100-tip boxes	712448	172.00
50-1250µL tips (orange/neutral)					
50-1250µL	Tip-Box 100-S	1000	10 filled 100-tip boxes	702450	182.00

Type of Tip

[A] 0.1-20µL nano-cap™ tips (dark gray/neutral)



[B] 0.5-20µL ultra-micro tips (dark gray/neutral)



[C] 2-200µL tips (yellow/neutral)



[D] 5-300µL tips (green/neutral)



[E] 50-1000µL tips (blue/neutral)




[F] 50-1250µL tips (orange/neutral)



NEW! Filter Pipette Tips

Applications that demand accuracy and precision as well as aerosol protection, such as PCR and isotope work, require the use of filter tips. BRAND filter tips combine the highest quality filters with premium tips to maximize tip performance by ensuring a great fit and a perfect seal. Tips are autoclavable at 121°C (250°F) and are made from dye-free resins.

- **Suitable for Critical Applications:** All BRAND filtered pipette tips are offered racked BIO-CERT® sterile—certified free of endotoxins, DNA, RNase and ATP—ensuring the purity and integrity of your sample. See page 19 for a description of BIO-CERT® purity.
- **Maintain Reaction Yield:** Tips are precision-molded from virgin polypropylene, with no added pigments, lubricants or trace-metals that can inhibit reactions. BRAND filters contain no sealing additives that can contaminate samples.
- **Ideal for Precise, Accurate Pipetting:** Tips feature clean, hydrophobic surfaces to reduce sample retention and increase accuracy and reproducibility. Tip orifice is precisely centered for directional accuracy.
- **Supplied in Easy-to-Use Packaging:** Available in a variety of packaging configurations to accommodate the requirements of any application.



BIO-CERT®

**Superior Quality Tips.
Aerosol Protection.**



BRAND filter pipette tips



BIO-CERT® Tip-Refill (sterile)

Filled tip inserts in low-waste, recyclable PET packaging. Designed to refill 60- or 96-tip format BRAND Tip-Boxes.



BIO-CERT® Filter Tip-Set (sterile)

Comes complete with one filled sterile Tip-Box S plus four BIO-CERT® filter tip refill inserts. 1000µL Tip-Set S comes with eight BIO-CERT® filter tip refill inserts.



BIO-CERT® Filter Tip-Box S & Filter Tip-Box 100-S (sterile)

Filled boxes of sterile tips. Tip-Boxes holding 96 racked filter tips may be refilled with BIO-CERT® Filter Tip-Refills. Tip-Box 100-S holds one hundred 1000µl tips in a disposable box.


For assistance in selecting the correct tips for your Transferpette® pipette, consult the ordering information section for your specific model. NOTE: filters reduce total tip capacity.

ORDERING INFORMATION


Capacity/ Volume	Packaging	# of tips	Description	Cat. No.	2011 List Price
NEW! 0.1-1µL nano-cap™ tips (dark gray/neutral)					
0.1-1µL	Filter Tip-Refill S	960	10 filled 96-tip inserts	702139	\$202.00
0.1-1µL	Filter Tip-Set S	480	5 filled 96-tip inserts, and one polypropylene box	702119	156.00
0.1-1µL	Filter Tip-Box S	960	10 filled polypropylene boxes with 96 tips each	702159	308.00
0.5-10µL tips (dark gray/neutral)					
0.5-10µL	Filter Tip-Refill S	960	10 filled 96-tip inserts	702140	181.00
0.5-10µL	Filter Tip-Set S	480	5 filled 96-tip inserts, and one polypropylene box	702120	103.00
0.5-10µL	Filter Tip-Box S	960	10 filled polypropylene boxes with 96 tips each	702160	224.00
NEW! 2-20µL tips (yellow/neutral)					
2-20µL	Filter Tip-Refill S	960	10 filled 96-tip inserts	712142	152.00
2-20µL	Filter Tip-Set S	480	5 filled 96-tip inserts, and one polypropylene box	712122	90.00
2-20µL	Filter Tip-Box S	960	10 filled polypropylene boxes with 96 tips each	712159	224.00
5-100µL tips (yellow/neutral)					
5-100µL	Filter Tip-Refill S	960	10 filled 96-tip inserts	702144	152.00
5-100µL	Filter Tip-Set S	480	5 filled 96-tip inserts, and one polypropylene box	702124	90.00
5-100µL	Filter Tip-Box S	960	10 filled polypropylene boxes with 96 tips each	702164	224.00
5-200µL tips (yellow/neutral)					
5-200µL	Filter Tip-Refill S	960	10 filled 96-tip inserts	702146	152.00
5-200µL	Filter Tip-Set S	480	5 filled 96-tip inserts and one polypropylene box	702126	90.00
5-200µL	Filter Tip-Box S	960	10 filled polypropylene boxes with 96 tips each	702166	224.00
50-1000µL tips (blue/neutral)					
50-1000µL	Filter Tip-Box 100-S	1000	10 filled 100-tip boxes	702168	286.00

Type of Tip

[I] 0.1-1µL nano-cap™ tips (Dark Gray/Neutral)



[J] 0.5-10µL tips (Dark Gray/Neutral)



[K] 2-20µL tips (Yellow/Neutral)



[L] 5-100µL tips (Yellow/Neutral)



[M] 5-200µL tips (Yellow/Neutral)



[N] 50-1000µL tips (Blue/Neutral)



NEW! Ultra Low Retention Tips

BRAND Ultra Low Retention tips feature a patented surface treatment to make the polypropylene of the pipette tips extremely liquid repellent thus reducing costly sample loss during pipetting.

- Ultra-hydrophobic tips are ideal for biological samples that contain detergents such as Triton® X-100, SDS, Tween etc.
- Special patented process does not use coatings that might lead to sample contamination
- Highly homogeneous surface is free of defect sites, preventing nonspecific binding



- High chemical resistance is ideal for working with solvents
- Autoclavable

The new ultra-hydrophobic Ultra Low Retention pipette and filter tips from BRAND.

Surface	Surface Tension
BRAND® PP Ultra Low Retention	9mN/m
PTFE	19mN/m
Silicone	21.5mN/m
Untreated PP	30mN/m
Water	72mN/m

A Persistent Problem

Limited volume samples, costly enzymes and reagents, and foaming liquids all need to be pipetted with high precision. However, close inspection of conventional pipette tips reveal that fluid droplets or a thin film of liquid remains in the pipette tip after dispensing. Significant factors leading to these residues are the different surface energy components and molecular composition of the pipette tips and the sample being pipetted.

The Dirty Secret of Detergents

In biological applications, such as PCR, PAGE or many other areas of protein chemistry, the sample must often be treated with detergent. An exact appraisal of the behavior of polypropylene surfaces in these applications is complicated. Detergents significantly decrease the surface tension of liquids, which hinders the optimal dispensing of the liquid from the pipette tip and diminishes pipette performance.

A Slick Solution

The new, ultra-hydrophobic BRAND pipette tips and filtered pipette tips solve the problem. The surfaces of the Ultra Low Retention tips are produced through a special, patented physiochemical process. The homogeneous, defect-free surface produced has extremely low surface tension. This significantly reduces sample loss and provides substantially higher reproducibility when working with critical media.



BRAND Ultra Low Retention pipette tips

ORDERING INFORMATION

Capacity/ Volume	Packaging	# of tips	Description	Cat. No.	2011 List Price
NEW! Ultra Low Retention Pipette Tips					
0.1-20µL nano-cap™ tips (dark gray/neutral)					
0.1-20µL	Tip-Box PP	480	5 polypropylene tip boxes filled with 96 tips each	703100	\$55.00
0.1-20µL	Tip-Box S	960	10 polypropylene tip boxes filled with 96 tips each, BIO-CERT® sterile	703122	153.00
0.5-20µL ultra-micro tips (dark gray/neutral)					
0.5-20µL	Tip-Box PP	480	5 polypropylene tip boxes filled with 96 tips each	703101	36.00
0.5-20µL	Tip-Box S	960	10 polypropylene tip boxes filled with 96 tips each, BIO-CERT® sterile	703123	100.00
0.5-20µL	Tip-Stack™	960	10 filled and stacked 96-tip inserts and a polypropylene box	703110	56.00
2-200µL tips (yellow/neutral)					
2-200µL	Tip-Box PP	480	5 polypropylene tip boxes filled with 96 tips each	703102	27.00
2-200µL	Tip-Box S	960	10 polypropylene tip boxes filled with 96 tips each, BIO-CERT® sterile	703124	95.00
2-200µL	Tip-Stack™	960	10 filled and stacked 96-tip inserts and a polypropylene box	703111	48.00
5-300µL tips (yellow/neutral)					
5-300µL	Tip-Box PP	480	5 polypropylene tip boxes filled with 96 tips each	703103	34.00
5-300µL	Tip-Box S	960	10 polypropylene tip boxes filled with 96 tips each, BIO-CERT® sterile	703125	95.00
50-1000µL tips (blue/neutral)					
50-1000µL	Tip-Box-100	500	5 polypropylene tip boxes filled with 100 tips each	703104	44.00
50-1000µL	Tip-Box S-100	1000	10 polypropylene tip boxes filled with 100 tips each, BIO-CERT® sterile	703126	115.00
50-1000µL	Tip-Stack™-100	1000	10 filled and stacked 100-tip inserts in two polypropylene boxes	703112	78.00
50-1250µL tips (orange/neutral)					
50-1250µL	Tip-Box-100	500	5 polypropylene tip boxes filled with 100 tips each	703105	56.00
50-1250µL	Tip-Box S-100	1000	10 polypropylene tip boxes filled with 100 tips each, BIO-CERT® sterile	703127	150.00
NEW! Ultra Low Retention Filter Pipette Tips					
0.1-1µL nano-cap™ filter tips (Dark Gray/Neutral)					
0.1-1µL	Filter Tip-Refill S	960	10 filled 96-tip inserts, BIO-CERT® sterile	703141	200.00
0.1-1µL	Filter Tip-Box S	960	10 polypropylene tip boxes filled with 96 tips each, BIO-CERT® sterile	703156	225.00
0.5-10µL filter tips (Dark Gray/Neutral)					
0.5-10µL	Filter Tip-Refill S	960	10 filled 96-tip inserts, BIO-CERT® sterile	703142	140.00
0.5-10µL	Filter Tip-Box S	960	10 polypropylene tip boxes filled with 96 tips each, BIO-CERT® sterile	703157	175.00
2-20µL filter tips (Yellow/Neutral)					
2-20µL	Filter Tip-Refill S	960	10 filled 96-tip inserts, BIO-CERT® sterile	703143	130.00
2-20µL	Filter Tip-Box S	960	10 polypropylene tip boxes filled with 96 tips each, BIO-CERT® sterile	703158	175.00
5-100µL filter tips (Yellow/Neutral)					
5-100µL	Filter Tip-Refill S	960	10 filled 96-tip inserts, BIO-CERT® sterile	703144	130.00
5-100µL	Filter Tip-Box S	960	10 polypropylene tip boxes filled with 96 tips each, BIO-CERT® sterile	703159	170.00
5-200µL filter tips (Yellow/Neutral)					
5-200µL	Filter Tip-Refill S	960	10 filled 96-tip inserts, BIO-CERT® sterile	703145	130.00
5-200µL	Filter Tip-Box S	960	10 polypropylene tip boxes filled with 96 tips each, BIO-CERT® sterile	703160	170.00
50-1000µL filter tips (Blue/Neutral)					
50-1000µL	Filter Tip-Box S-100	1000	10 polypropylene tip boxes filled with 100 tips each, BIO-CERT® sterile	703161	220.00

Type of Tip

[A] 0.1-20µL nano-cap™ tips (dark gray/neutral)



[B] 0.5-20µL ultra-micro tips (dark gray/neutral)



[C] 2-200µL tips (yellow/neutral)



[D] 5-300µL tips (green/neutral)



[E] 50-1000µL tips (blue/neutral)




[F] 50-1250µL tips (orange/neutral)



Type of Tip

[I] 0.1-1µL nano-cap™ tips (Dark Gray/Neutral)



[J] 0.5-10µL tips (Dark Gray/Neutral)



[K] 2-20µL tips (Yellow/Neutral)



[L] 5-100µL tips (Yellow/Neutral)



[M] 5-200µL tips (Yellow/Neutral)



[N] 50-1000µL tips (Blue/Neutral)

