



Lab Vacuum "VACUUM PUMPS NOW!" PROGRAM FOR 2012

Our most popular vacuum pumps can ship within one business day. Ask for "direct-shipment."

Buy products from leading lab dealers between 1/1/12 and 12/31/12.
For product information, visit www.brandtech.com, or call Customer Service at 888-522-2726.



Your vacuum pump just died, and with it your vacuum oven, freeze dryer, rotary evaporator, concentrator or gel dryer. You need a new pump, now. Where do you turn?

BrandTech Scientific announces its "Vacuum Pumps...Now!" program. We know you need your pump right away, so we will guarantee* that your order for any one of the pumps below will be shipped within one business day of your order. We'll even do our best to overnight the shipment* if your project is that critical.

So what do you get besides speedy delivery?

Select one of our oil-free pumps with a fluoropolymer, corrosion-resistant flowpath, and you'll get a whisper-quiet pump that will never need an oil change and typically needs its first service only after 10,000-15,000 hours of operation. That's years in most applications! And these oil-free pumps will handle nearly any lab evaporative application except freeze drying.

Choose one of our rotary vane pumps, and you'll get the very latest in oil-sealed pump technology. These totally redesigned pumps have a very small footprint, extremely low noise and vibration, great corrosion resistance, and vacuum pumping capacity of some competitive pumps twice their size.

Vacuum pumps for every application

BrandTech offers a complete line of vacuum pumps for the laboratory, from free-standing pumps to our unique self-regulating VARIO™ pumping systems. Many of these pumps are in stock for immediate delivery, and can be purchased through leading laboratory supply dealers. If you would like help choosing the best replacement pump for your application, call our factory-trained customer service department, or try the interactive Vacuum Pump Selection Guide on our website.

Be up and running in no time, with a quiet VACUUBRAND® pump from BrandTech!

* Request Vacuum Pumps NOW service and we guarantee that BrandTech will ship the above pump models direct from our warehouse within one business day or we will award a customer credit of \$150 worth of lab products direct from BrandTech. Guarantee applies to single-pump orders only direct-shipped from BrandTech. Other models may be available for immediate shipment. Customer pays shipping. Offer may be discontinued without notice. Visit www.brandtech.com for offer details and promotional specials.

BrandTech Scientific
Your ACE in the RAT RACE!

Product appearance, catalog numbers, prices, specifications, and technical information are subject to change without notice. Offer only valid to customers in the United States of America. Not valid for resellers. Void where prohibited by law or company policy. This offer may not be combined with any other offer.

Playing cards © 2008 BrandTech Scientific, Inc.

VACUUBRAND® vacuum pumps...NOW!

Model	Ultimate Vacuum, Torr	Flowrate, cfm (lpm)	BrandTech Cat. No.
Oil-free Chemistry-design diaphragm vacuum pumps			
ME1	75	0.5 (14)	721003
ME1C	75	0.5 (14)	721103
MZ2C NT	5	1.4 (38)	732303
MD1C	1.5	0.82 (23)	696613
Two-stage rotary vane vacuum pumps			
RZ2.5	1.5×10^{-3}	1.6 (47)	698123
RZ6	1.5×10^{-3}	4 (133)	698133