



BrandTech® Scientific

# PROMOTIONS FOR 2012

BRANDTECH®  
SCIENTIFIC, INC.

Toll Free 888-522-2726  
www.brandtech.com

Promotion redemption forms available at [www.brandtech.com](http://www.brandtech.com)

### High Card

- **HandyStep® electronic** — Buy a *HandyStep® electronic* repeating pipette, get 3 boxes of PD-Tip™ syringe tips (non-sterile or BIO-CERT® pure) FREE!

### One Pair

- **Transferpette® S pipettes** — Buy 2 *Transferpette® S* single channel pipettes, get one FREE!
- **Transferpette® multichannel pipettes** — Buy 2 *Transferpette®* manual multichannel mechanical pipettes (classic or *S* style), get one FREE!
- **Transferpette® electronic pipettes** — Buy 2 *Transferpette® electronic* pipettes (single or multichannel) get one single channel *Transferpette® electronic* FREE!
- **UV-Cuvettes** — Buy 2 packages of UV-Cuvettes, get a package of PCR tubestrips and caps FREE!

### Two Pair

- **macro™** pipette controller — Buy 2 macros, get two macros FREE!

### Three of a Kind

- **accu-jet® pro** pipette controller — Buy 3 *accu-jet® pros*, get one FREE!
- **Dispensette® dispensers** — Buy 3 *Dispensette®* dispensers, get one FREE!

### Full House

- **HandyStep® S & PD-Tip™ syringe tips** — Buy \$400 worth of PD-Tip™ syringe tips and get a *HandyStep® S* repeating pipette FREE!

### All-In!

- **Any BrandTech product** — Buy \$500 worth of BrandTech products, get \$100 (list price) in BrandTech products FREE!

### Wild Card

- **FREE Lab Rat T-Shirt** — With any BrandTech purchase over \$200!
- **VACUUBRAND® vacuum pump** — Buy a vacuum pump or system, get \$500\* (list price) in BrandTech products for FREE! (\*\$250 with purchase of Tyro12™, \*ME1 & ME1C excluded).

### Call our Bluff!

- **“Vacuum Pumps NOW!”** — Buy one of our 5 top-selling vacuum pumps, get it shipped within one business day!



**BrandTech® Scientific — Your ACE in the RAT RACE!**

Offers only valid to customers in the United States of America. Not valid for resellers. Void where prohibited by law or company policy. May not be combined with other offers.

Playing cards © 2008 BrandTech Scientific, Inc.



**Bottletop Dispensers & Burettes**



**Manual, Electronic & Repeating Pipettes**



**Pipette Controllers**



**Laboratory Vacuum**



**Life Science Plastics**




**Volumetric & General Lab Plastics**



**Stainless Support Jacks**

## BRAND<sup>®</sup> plates<sup>®</sup> microplates for cell culture, immunoassay & general use

- 96-, 384- and 1536-well
- Non-treated & treated surfaces for immunoassay & cell culture
- Transparent, white, black, clear bottom and UV-transparent
- C, F, V and U-shaped wells
- Deep Well Plates
- Barcoding on request
- Physical treatments that eliminate refrigeration and shelf-life concerns
- Samples available



## BRAND PLT<sup>™</sup> Pipette Leak Tester




- Validate pipette function in seconds
- Use routinely between calibrations
- Excellent diagnostic tool

The BRAND<sup>®</sup> PLT<sup>™</sup> Pipette Leak Tester is an entirely new class of device intended to verify the performance of pipettes between regularly scheduled calibration. It quantifies instrument leakage—the cause of 80% of instrument failures, and provides a pass/fail result based on whether it will provide a gravimetrically detectable error. By testing with or without tips, statically or dynamically, the PLT<sup>™</sup> helps pinpoint the source of an instrument failure quickly and effectively.



## BRAND HandyStep<sup>®</sup> electronic Repeating Pipette

- Motorized, push-button operation eliminates thumb fatigue
- Choose virtually any dispensing volume between 1µL and 50mL
- Instrument “learns” operator’s dispensing rhythm without programming
- Compatible with most leading syringe tips, including dealer private label tips
- Automatically recognizes BRAND PD-Tip<sup>™</sup> syringe tips upon insertion
- **Mechanical HandyStep<sup>®</sup> S also available!**



## BRAND PD-Tip<sup>™</sup> Syringe Tips / Repeating Pipette Tips

- Manufactured from high-quality PP, HDPE and LCP plastics
- Encoded to be automatically recognized by the HandyStep<sup>®</sup> electronic and the Rainin<sup>®</sup> AutoRep<sup>™</sup> E repeating pipettes. Compatible with the BRAND HandyStep<sup>®</sup> S, Eppendorf<sup>®</sup> Repeater<sup>®</sup>, and other repeating pipettes using standard-sized tips

## VACUUBRAND<sup>®</sup> ME1 and ME1C Vacuum Pumps Compact vacuum pumps for filtration and SPE

The new VACUUBRAND<sup>®</sup> ME1 and ME1C oil-free diaphragm pumps offer a compact, high performance and easy to-use solution which is perfect for single and multiple filtrations as well as solid phase extraction.



- Compact small footprint conserves valuable bench space
- Whisper Quiet minimizes background noise in the lab
- Long Diaphragm Life reduces maintenance-related downtime
- Top Mounted Power Switch simplifies operation
- High performance 141lpm pumping speed and 100mbar ultimate vacuum

Both models feature robust PTFE diaphragms and valves for enhanced chemical resistance. The ME1C features a full fluoropolymer flowpath to handle more aggressive vapors. An optional vent control valve with dial gauge enables variable fine adjustment of the pumping speed.

## VITLAB<sup>®</sup> PFA, Volumetric and General Lab Plastics



PFA Labware for Trace Analysis



Erlenmeyer Flasks



Light Shielding



Class A & B Graduated Cylinders and Class A & B Volumetric Flasks with stoppers or screw caps

Details at [www.brandtech.com](http://www.brandtech.com) or call for a catalog of our line of PFA, Volumetric & General Laboratory plastic products, including bottles, beakers, sample vials, funnels and scoops.



Beakers, Pitchers, Scoops, Funnels and more!

Custom imprinting available

VITLAB

**BRANDTECH<sup>®</sup>**  
SCIENTIFIC, INC.

Toll Free 888-522-2726  
[www.brandtech.com](http://www.brandtech.com)



Lab Rats Trust BrandTech<sup>®</sup>!