

# Reliable Tools for Industrial Labs

## The Standard in Bottletop Dispensers

### Dispensette® S dispensers

- Rugged: no seals to wear out!
- Suitable for most laboratory reagents
- Simple to clean & calibrate
- Recirculation valve reduces splashing and helps conserve reagents
- Autoclavable to 121 °C
- Adapters for common bottle threads included
- Optional flexible discharge tube simplifies serial dispensing



Dispensette S	recirculation valve	Recirculation valve
Volume, mL	Cat. No.	Cat. No.
<b>Dispensette S Digital</b>		
0.1-1	4600310	4600311
0.2-2	4600320	4600321
0.5-5	4600330	4600331
1-10	4600340	4600341
2.5-25	4600350	4600351
5-50	4600360	4600361
<b>Dispensette S, Analog-adjustable</b>		
0.1-1	4600100	4600101
0.2-2	4600120	4600121
0.5-5	4600130	4600131
1-10	4600140	4600141
2.5-25	4600150	4600151
5-50	4600160	4600161
10-100	4600170	4600171

Fixed volume models available. For strong acids and nonpolar organics, ask for the BRAND® Dispensette® S Organic dispensers.

## Faster, Easier Titrations with Class A Precision

### Titrette® bottletop burette



- Meets Class A standards for glass burettes!
- Popular for sulfur dioxide (SO<sub>2</sub>) analyses, acid-base titrations, solvent analyses and more
- Suitable for NaOH, KOH, KI, H<sub>2</sub>SO<sub>4</sub>, KMnO<sub>4</sub>, Iodine and more
- Eliminates meniscus reading errors
- Models available with direct data transmission to PC

Description	Cat. No.
10 mL (standard)	4760141
25 mL (standard)	4760151
50 mL (standard)	4760161

## Convenient & Economical

### Disposable PS and PMMA Cuvettes


- Consistent reliable results—manufactured from the highest quality virgin resins
- Accommodate most applications - available macro and semi-micro sizes
- Easy to use—no more cleaning



## Class A, USP, Certified Volumetric Glassware

Precision analysis requires measuring instruments with consistently high precision. BLAUBRAND® volumetric instruments represent the highest degree of technical perfection. They are made of glass which is highly resistant to nearly all chemical substances.

- Volumetric instruments are individually calibrated (USP batch certificate supplied)
- High quality inks and annealing process provide durable graduations and long service life
- Thermal stress in the glass blanks is eliminated prior to adjustment so there are no permanent volume changes after heating, up to 250 °C



### BLAUBRAND® graduated cylinder, tall form, Class A, USP, Certified, Boro 3.3

Volume (mL)	Height (mm)	Qty./Pack	Cat. No.
5	115	2	32805
10	140	2	32808
25	170	2	32820
50	200	2	32828
100	260	2	32838
250	335	2	32848
500	365	2	32854
1000	465	1	32862
2000	505	1	32864

# Reliable Tools Industrial Labs

## Precise, Slim, Light & Dependable Repeater

### HandyStep® S Repeating Pipette



The HandyStep S provides unprecedented versatility in a mechanical repeating pipette. It offers 59 different volume settings with accuracy, precision and dependability.

- Comfortable, slim design is lightweight, well-balanced and offers easy one-handed volume selection
- Flexible; accepts third party tips – no more being locked into a sole source
- Tip ejection system allows tip removal without user handling
- HandyStep touch and HandyStep touch S models also available

Description	Cat. No.	
HandyStep S Repeating Pipette	705110	
HandyStep touch Repeating Pipette	705200	
<b>PD-Tip™ // size</b>	<b>Non-sterile tips</b>	<b>BIO-CERT® Pure*</b>
0.1 mL (100 tips)	705700	705730
0.5 mL (100 tips)	705702	705732
1.0 mL (100 tips)	705704	705734
1.25 mL (100 tips)	705706	705736
2.5 mL (100 tips)	705708	705738
5 mL (100 tips)	705710	705740
10.0 mL (100 tips)	705712	705742
12.5 mL (100 tips)	705714	705744
25 mL (50 tips)	705716	—
25 mL (25 tips)	—	705746
50 mL (25 tips)	705718	705748
Multi-pack (20 tips each 0.5, 1, 1.25, 2.5, 5, 10 & 12.5 mL)	705720	—

\*BIO-CERT® is certified sterile, and DNA-, RNase-, ATP- and endotoxin-free.

## Lightweight and Rugged

### Transferpette® S Adjustable Volume Single Channel Pipettes



Capacity, µL	Cat. No.
0.1-1	705868
0.1-2.5	705869
0.5-10	705870
2-20	705872
5-50	705873
10-100	705874
20-200	705878
100-1000	705880
500-5000 (0.5-5 mL)	705882
1000-10000 (1-10 mL)	705884

### Transferpette® S -8/-12 Multichannel Pipettes



Capacity, µL	8-Channel	12-Channel
0.5-10	705900	705920
5-50	705906	705926
10-100	705908	705928
20-200	705910	705930
30-300	705912	705932

Transferpette S -8/-12 pipettes supplied with performance certificate, operating manual, shelf mount, two tip boxes, reagent reservoir, pack of replacement O-rings. One-year warranty.

## Comfortable, Convenient Meniscus Control

### macro™ Pipette Controller



- Replaces old-fashioned pipette bulbs with simple, one-finger control
- Rapid or sensitive filling; fast or drop-wise dispensing
- Fits most pipettes 0.1-200 mL

Description	Cat. No.
Gray	26200
Blue	26202
Magenta	26203
Green	26201