

Discharge Tubes for seripettor®; and seripettor®; pro



Price: \$39.60 - \$81.80

Product Categories: [Bottletop Dispensers Parts & Accessories](#)

Product Page:

<https://www.brandtech.com/product/discharge-tubes-for-seripettor-and-seripettor-pro/>

Product Variants

- Discharge Tubes for seripettor® and seripettor® pro - seripettor® pro, 10mL ()
- Discharge Tubes for seripettor® and seripettor® pro - seripettor®, 2mL ()
- Discharge Tubes for seripettor® and seripettor® pro - seripettor®, 10, 25mL ()
- Discharge Tubes for seripettor® and seripettor® pro - seripettor® pro, 2mL

()

- Discharge Tubes for seripettor® and seripettor® pro - seripettor® pro, 25mL

()

- Discharge Tubes for seripettor® and seripettor® pro - ADAPTER for seripettor® pro discharge tube, 2, 10, 25mL ()

Product Description

Discharge tubes for the seripettor® pro have an integrated discharge valve and require an adapter for attachment to the valve body.

[SafetyPrime Valves](#) for the seripettor® pro and [valve sets](#) for the seripettor® are sold separately.

Product Attributes

- Dimensions: N/A

- product: ADAPTER for seripettor®; pro discharge tube, Dispensette®; III, Dispensette®; Organic, Dispensette®; TA, Pt-Ir, Dispensette®; TA, Ta, seripettor®, seripettor® pro

- volume: 0.5, 1, 2, 5, 10mL, 0.5, 1, 2, 5mL, 10, 25mL, 10mL, 2, 10, 25mL, 25, 50, 100mL, 25mL, 2mL, 5, 10mL