

Griffin Beakers



\$13.00 - \$207.00][15, 117, 84][14]

Product Categories: [Plasticware](#)

Product Page: <https://www.brandtech.com/product/griffin-beakers/>

Product Variants

- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 5000mL
- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 10mL
- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 25mL
- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 50mL

- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 100mL

- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 150mL

- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 250mL

- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 400mL

- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 500mL

- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 600mL

- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 800mL

- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 1000mL

- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 2000mL

- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 3000mL

- Griffin Beakers ()material: PP, options: Molded graduations, volume: 25mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 50mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 100mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 150mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 250mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 400mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 500mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 600mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 1000mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 2000mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 3000mL
- Griffin Beakers ()material: PP, options: Molded graduations, volume: 5000mL
- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening, volume: 10mL
- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening, volume: 25mL
- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,

volume: 50mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,
volume: 100mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,
volume: 150mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,
volume: 250mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,
volume: 400mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,
volume: 500mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,
volume: 600mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,
volume: 800mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,
volume: 1000mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,
volume: 2000mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening,
volume: 3000mL

- Griffin Beakers ()material: PP, options: Molded graduations w/blue screening, volume: 5000mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 25mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 50mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 100mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 150mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 250mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 400mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 500mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 600mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 1000mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 2000mL
- Griffin Beakers ()material: PMP, options: Molded graduations, volume: 3000mL
- Griffin Beakers ()material: PMP, options: Printed red scale, not raised, volume: 5000mL

- Griffin Beakers ()material: ETFE, options: Printed black scale, not raised,
volume: 25mL

- Griffin Beakers ()material: ETFE, options: Printed black scale, not raised,
volume: 50mL

- Griffin Beakers ()material: ETFE, options: Printed black scale, not raised,
volume: 100mL

- Griffin Beakers ()material: ETFE, options: Printed black scale, not raised,
volume: 250mL

- Griffin Beakers ()material: ETFE, options: Printed black scale, not raised,
volume: 400mL

- Griffin Beakers ()material: ETFE, options: Printed black scale, not raised,
volume: 500mL

- Griffin Beakers ()material: ETFE, options: Printed black scale, not raised,
volume: 600mL

- Griffin Beakers ()material: ETFE, options: Printed black scale, not raised,
volume: 1000mL

Product Summary

Product and ordering information for Griffin Beakers. Available in a variety of plastics to fit your application: ETFE, PMP or PP. Meets the requirements of ISO 7056.

Product Description

These Griffin beakers are available in a variety of plastics to fit your application. ETFE (ethylene tetrafluoroethylene) provides high chemical resistance and very high thermal stability (-100 to +150°C). PMP (polymethylpentene) provides glass-like clarity, while PP (polypropylene) is highly transparent and well-suited for general purpose use. Meets the requirements of ISO 7056 (except V110204-V111004). To preserve markings on products with colored overprinting (cat. no. V605081-V617081, V60503-V61703, V110204-V111004), cleaning at no higher than 60°C is recommended. For autoclaving, the design with raised graduations (cat. no. V606941-V617941, V60695-V61795) is recommended.

Product Attributes

- Dimensions: N/A
- material: ETFE, PMP, PP
- options: Molded graduations, Molded graduations w/blue screening, Printed black scale, not raised, Printed red scale, not raised
- volume: 1000mL, 100mL, 10mL, 150mL, 2000mL, 250mL, 25mL, 3000mL, 400mL, 5000mL, 500mL, 50mL, 600mL, 800mL

Product Gallery



