

Mucosal Universal Detergent



Price: \$166.60

SKU: 44003B

Product Categories: [Labware](#)

Product Page: <https://www.brandtech.com/product/mucosal-universal-detergent/>

Product Summary

Product and ordering information for Mucosal® Universal Detergent. Universal cleaner for laboratory equipment and instruments - glass, porcelain, plastic, rubber and metal. Non-toxic, non-aggressive, and biodegradable.

Product Description

A liquid, alkaline, high-performance concentrate for immersion and ultrasonic baths, Mucosal® is a universal cleaner for laboratory equipment and instruments made of glass, porcelain, plastic, rubber and metal. Items are cleaned

in ~10-30 minutes, without mechanical help. Mucaso[®] replaces chromosulfuric acid, is free from chlorine and caustic alkalis, is non-toxic, non-aggressive, and biodegradable.

Ingredients:

> 30% phosphates

< 5% anionic surfactants, amphoteric surfactants, complexing agents, corrosion inhibitors, auxiliary agents

Product Gallery

