

NEW! Fluorimetry Cuvettes



Price: \$20.40 - \$408.60

Product Categories: [Cuvettes](#), [Cuvettes](#)

Product Page: <https://www.brandtech.com/product/new-fluorimetry-cuvettes/>

Product Variants

- NEW! Fluorimetry Cuvettes - PS, Pack of 100 ()
- NEW! Fluorimetry Cuvettes - PS, Pack of 500 ()
- NEW! Fluorimetry Cuvettes - UV, Pack of 100 ()
- NEW! Fluorimetry Cuvettes - UV, Pack of 500 ()

Product Summary

Product and ordering information for Fluorimetry cuvettes. Four-sided clear disposable cuvettes for fluorescence applications; Polystyrene for 340nm-800nm or proprietary polymer for 230nm-800nm.

Product Description

NEW! Disposable macro cuvettes with 4 optically clear sides from BRAND suitable for use in fluorescence applications. Available in high quality polystyrene for measurements from 340nm to 800nm, as well as a proprietary polymer with ultra-low autofluorescence for precise measurements in both the visible and UV-range (230-800nm). Cuvettes have a standard 10mm light path for use in all standard fluorimeters. Packaging is 100 per tray in trays made of low-dust,

non-scratching polystyrene, and are grouped by manufacturing mold cavity to ensure the lowest variation in extinction coefficient.

Product Attributes

- Dimensions: N/A
- material: PS, UV
- options: Pack of 100, Pack of 500