

QuikSip™ Bottletop Aspirator



Price: \$334.60

SKU: 4723150

Product Categories: [Fluid Aspiration](#), [Fluid Aspiration](#)

Product Page: <https://www.brandtech.com/product/quiksip-bottletop-aspirator/>

Product Summary

Product and ordering information for the QuikSip Bottletop Aspirator: Convenient bottletop aspirator for removing supernatants, cell culture media and other biological liquids.

Product Description

The QuikSip™ is a convenient bottletop aspirator that works without any external vacuum source, allowing it to be easily transferred between workstations. Each plunger stroke powers up to 25mL of aspiration, and its unique

thumb-operated controller is compatible with disposable pipette tips, micropipettes, and glass Pasteur pipettes.

□ Product Features:

No vacuum required.

Unique thumb-operated controller compatible with disposable pipette tips, micropipettes and glass Pasteur pipettes.

Autoclavable for sterile applications.

Mounts on most common reagent bottles (45mm thread; 45/33mm and 45/S40mm adapters supplied).

□

Product Attributes

- Dimensions: 12.75 × 8.125 × 3 in
- Weight: 1.5 lbs