

## seripettor® Bottletop Dispenser



Price: \$132.00 - \$167.00

**Product Categories:** [Bottletop Dispensers](#)

**Product Page:** <https://www.brandtech.com/product/seripettor/>

### Product Variants

- seripettor® Bottletop Dispenser - 0.2-2mL ( )
- seripettor® Bottletop Dispenser - 1-10mL ( )
- seripettor® Bottletop Dispenser - 2.5-25mL ( )

### Product Summary

Product and ordering information for the seripettor wiping-seal bottletop dispenser

suitable for many reagents commonly used in biological and chemical laboratories.

## Product Description

The seripettor® is an economical wiping seal dispenser suitable for many commonly used reagents in biological, clinical, medical, and chemical laboratories. Buffers, bases, culture media, even sodium azide solutions, can all be dispensed easily and precisely. Sterile dispensing of media or agar can be performed with additional accessories.

It simplifies dispensing with its notched volume selector for quick and exact volume adjustment. It also features a spring-driven refill function for semi-automatic operation, and the innovative design of the seripettor permits simple disassembly and rapid replacement of wearing parts, extending instrument life.

## Product Attributes

- Dimensions: N/A
- volume: 0.2-2mL, 1-10mL, 2.5-25mL