



VITsafe™ Safety Wash Bottles, LDPE, spray tube PP

- Durable safety imprint
- Material name in German, English, French, and Spanish
- Chemical formula, CAS No. hazard pictogram, signal word, risk phrases (H phrases), safety phrases (P phrases), as well as NFPA code

Reagent	Cap Color	Volume mL	Diameter mm	Height* mm	Qty./Pack	Cat. No.	2020 List Price
Distilled Water	Blue	500	76	166	12	V1352819	\$209.00
Ethanol	Orange	500	76	166	12	V1352869	209.00
Acetone	Red	500	76	166	12	V1452829	209.00
Heptane	Red	500	76	166	6	V1352899	106.40
Hexane	Red	500	76	166	6	V1352909	106.40
Xylene	Red	500	76	166	6	V1352959	106.40

*Height without spray tube



Dropping Bottles, LDPE, with dropper screw cap, HDPE



Volume mL	Thread GL	Diameter mm	Height mm	Qty./Pack	Cat. No.	2020 List Price
50	18	37	129	24	V94587	\$68.20
100	18	43	155	24	V94687	76.60
250	25	58	183	12	V94987	56.60
500	25	74	228	12	V95087	74.40
1000	32	92	269	12	V95187	100.80

Reagent Bottles, PP, with NS stoppers, PP

- Good chemical resistance, ideal for long-term storage of liquids
- Autoclavable at 121° C according to DIN EN 285



Volume mL	Stopper* NS	Height mm	Diameter mm	Qty./Pack	Cat. No.	2020 List Price
100	14/23	106	52	20	V100394	\$72.40
1000	29/32	213	108	10	V100694	96.60

*Stopper type A: With square knob and red core.

Reagent Bottles, PP, with screw caps, PP



- Screw caps help protect against outside contamination and leakage
- Autoclavable to 121° C according to DIN EN 285, except for V100989 and V100889



Volume mL	Thread GL	Height mm	Diameter mm	Qty./Pack	Cat. No.	2020 List Price
100	18	100	52	20	V100389	\$64.00
250	25	132	70	20	V100489	85.00
500	25	165	87	10	V100589	61.00
1000	32	202	108	10	V100689	80.80
2000	32	245	131	6	V100789	100.80
5000	45	315	178	1	V100889*	29.40
10000	63	394	222	1	V100989*	73.40

*Bottle V100889 includes 1 handle; bottle V100989 includes 2 handles.

VITLAB® General Laboratory Products

Wide Mouth Reagent Bottles, PP, with screw caps, PP

- Screw caps help protect against outside contamination and leakage
- Wider mouth for easier filling
- Autoclavable at 121° C according to DIN EN 285



Volume mL	Thread GL	Height mm	Diameter mm	Qty./Pack	Cat. No.	2020 List Price
100	32	96	55	20	V101589	\$85.00
250	45	132	73	20	V101689	108.00
500	45	172	87	10	V101789	69.20
1000	63	204	108	10	V101889	125.00
2000	63	243	131	6	V101989	94.40



Bottles, LDPE, with screw caps, PP

- Screw caps help protect against outside contamination and leakage



50	18	85	37	24	V94589	65.00
100	18	144	43	24	V94689	69.20
250	25	135	65	12	V94989	48.40
500	25	180	74	12	V95089	69.20



Wide Mouth Bottles LDPE, with screw caps, PP

- Screw caps help protect against outside contamination and leakage
- Wider mouth for easier filling



250	45	146	58	12	V93789	55.60
500	45	166	76	12	V93989	76.60
1000	63	226	91	12	V94189	109.00



Square Bottles, HDPE, with screw caps, LDPE

- Screw caps help protect against outside contamination and leakage
- Space-saving smaller footprint with square-shaped bottle

Volume mL	Thread GL	Height mm	Size mm	Qty./Pack	Cat. No.	2020 List Price
100	25	76	43 x 43	24	V91789	\$77.60
250	28	80	80 x 80	24	V91989	104.00
500	32	106	90 x 90	12	V92089	69.20
1000	32	187	80 x 80	12	V92189	94.40




Wide Mouth Square Bottles, HDPE, with screw caps, LDPE

- Screw caps help protect against outside contamination and leakage
- Space-saving smaller footprint with square-shaped bottle

100	32	78	46 x 46	24	V92489	102.80
250	50	83	80 x 80	24	V92689	116.40
500	65	120	90 x 90	12	V92789	87.00
1000	65	168	90 x 90	12	V92889	109.00



VITLAB® General Laboratory Products



Forceps

Length mm	Qty./Pack	Cat. No.	2020 List Price
--------------	-----------	----------	--------------------

Forceps, PMP with pointed ends

115	10	V67895	\$28.40
145	10	V67995	33.60



Forceps, POM with rounded ends

115	5	V68099	13.20
145	5	V68199	16.00
180	5	V68299	30.00
250	5	V68399	47.40



Scoops

Volume mL	Length mm	Qty./Pack	Cat. No.	2020 List Price
--------------	--------------	-----------	----------	--------------------

Measuring Scoops, PP, volume marked on handle


2	60	12	V39194	\$6.60
5	82	12	V39294	8.00
10	100	12	V39394	10.40
25	135	12	V39494	11.00
50	160	12	V39594	21.00
100*	200	12	V39694	27.40
250	260	6	V39794	34.60
500	315	6	V39894	48.40
1000	385	6	V39994	59.80

*Colors available – inquire



Measuring Scoops, HDPE

15	115	12	V40093	11.60
25	135	12	V40193	17.00
65	185	12	V40293	22.00
110	215	12	V40393	31.40
150	250	12	V40493	36.80
350	310	6	V40593	35.60
750	350	6	V40693	58.80
1250	400	6	V40793	82.00



Slide Staining Kit, POM

Description	Size mm	Height mm	Qty./Pack	Cat. No.	2020 List Price
Universal Slide Staining Kit*	100 x 87	51	5	471800	\$193.00
Staining jar	100 x 87	51	5	471500	158.00
Staining rack	91 x 79	38	5	471400	66.20

*Universal Kit 471800 includes one staining jar and a staining rack for 25 slides, 76 x 26 mm.

VITLAB 

Please see page 130 for description of various plastics.

VITLAB® General Laboratory Products

Funnels, PP

- Handle with loop for hanging

Diameter mm	Approx vol. mL	Length mm	Stem Diameter mm	Qty./pack	Cat. No.	2020 List Price
30	5	45	1.5	24	V40894	\$21.00
30	6	47	4	24	V41094	21.00
40	14	65	4	24	V41194	25.20
50	32	85	7	24	V41294	31.40
75	88	108	7	12	V41394	20.00
100	222	155	8	12	V41494	37.80
120	342	180	11	12	V41594	47.40
150	817	220	15	12	V41694	63.00



Funnels

Large Funnels, PP or HDPE*

- Integrated handle (except V41994)

200 (PP)	1300	200	22	6	V41794	72.40
250 (PP)	3200	260	30	6	V41894	126.00
350 (PP)	12500	440	35	1	V41994	44.00
400 (HDPE)*	12500	365	42	1	V42294	112.40
430 (HDPE)*	17500	420	37	1	V42393	149.00

*HDPE Large Funnels V42294 and V42393 are not food-safe products.



Large funnels

Powder Funnels, PP

65	70	15	10	V70794	27.00
80	75	21	10	V70894	38.80
100	92	24	10	V70994	41.00
120	105	27	10	V71094	50.40
150	138	28	5	V71194	32.60



Powder funnels

Sample Vials, PP, with snap caps, LDPE

Volume mL	Height mm	Diameter mm	Qty./pack	Cat. No.	2020 List Price
5	25	20	25	V68594	\$48.40
18	57	22	25	V68894	62.00
50	97	30	10	V69194	32.20
160	110	50	10	V69294	37.60



Sample Vials, with snap caps

Sample Vials with screw caps

Sample Vials, PP, with screw caps, PP

Volume mL	Height mm	Diameter mm	Thread GL	Qty./pack	Cat. No.	2020 List Price
60	90	38	40	10	V130394	\$111.40
180	112	54	56	10	V130594	162.80



Please see page 130 for description of various plastics.



Bucket, PP



Bucket, HDPE

Buckets, PP, with handles, pouring lips, and graduations.



- With handle and lip for ease of pouring
- Markings at 1L graduations

Volume L	Height mm	Diameter mm	Qty./Pack	Cat. No.	2020 List Price
12	330	310	1	V96694	\$49.40
15	370	310	1	V96794	55.60

Buckets, HDPE, with handles



- With carrying handle
- Markings at 1L graduations

5	240	250	1	V96093	24.40
10	300	290	1	V96393	29.40
Lid for 5L bucket (LDPE)	—	—	1	V96293	6.60
Lid for 10L bucket (LDPE)	—	—	1	V96593	7.60



Imhoff Sedimentation Cones, SAN

- Transparent
- Graduated to 1000 mL
- Screw closure allows drainage

1000 mL Cone	1	388000	36.20
Rack, PMMA, (for 2 cones) 150x300x290 mm	1	388050	311.40

VITLAB® custom imprinting

VITLAB® has specialized capabilities in high-precision printing on plastic products with superior chemical and break resistance. Many VITLAB® products can be imprinted with logos or custom graduations for OEM or promotional use in relatively small quantities. Custom imprinted VITLAB® products are ideal for handling granulates, powders or liquids commonly used in agriculture, laboratories, and in the food, beverage and chemical industries. From multi-color logos to specialized measuring scales, VITLAB® can imprint plastic labware to meet your individual requirements and specifications. Custom imprinted labware by VITLAB® is of the highest quality, and offers a unique and functional means of presenting your company brand. Beakers, pitchers and scoops can also be made in colored plastic.

Contact BrandTech® Scientific for more information.

VITLAB 

