

Description	Cat. No.	2020 List Price
<b>BVC Accessories</b>		
4 liter collection bottle, polypropylene, with aspiration cap	<a href="#">20635810</a>	\$478.60
2 liter collection bottle, glass, with aspiration cap	<a href="#">20635809</a>	537.60
2 liter bottle adapter for BVC control and BVC professional (included in BVC control G and BVC professional G)	<a href="#">20635839</a>	83.40
VHC <sup>PRO</sup> Vacuum Hand Controller with 2.5 meters of tubing	<a href="#">20688061</a>	257.20
8-channel pipette tip adapter, with ejector for VHC <sup>PRO</sup>	<a href="#">20635679</a>	455.60
Pipette tip adapter w/ ejector for VHC <sup>PRO</sup>	<a href="#">20635638</a>	83.40
Stand for VHC <sup>PRO</sup>	<a href="#">20635680</a>	158.40
Connector, to add second VHC <sup>PRO</sup> to BVC bottle aspiration cap	<a href="#">20699943</a>	67.20
Aerosol filter, 0.2µm	<a href="#">20638266</a>	52.20
Quick coupling to connect VHC <sup>PRO</sup> to BVC bottle aspiration cap	<a href="#">20635807</a>	178.40
Quick coupling to connect bottle to pump	<a href="#">20635808</a>	156.40
BVC Shuttle - allows easy movement of BVC under Biosafety cabinet	<a href="#">20696880</a>	222.20



VHC<sup>PRO</sup> with stand



VHC<sup>PRO</sup> 8-channel adapter



BVC Shuttle

## QuikSip™ Bottletop Aspirator

The QuikSip™ is a convenient bottletop aspirator that removes supernatants, cell culture media, and other biological liquids from dishes, flasks, culture bottles, centrifugation tubes, multiwell plates, and more. It is small, portable, and does not require an external vacuum or power source. The QuikSip™ includes a unique thumb-operated controller that is compatible with disposable pipette tips, micropipettes, and glass Pasteur pipettes.

### Operating limits and exclusions

Liquids unsuitable for use with the QuikSip™ bottletop aspirator include, but are not limited to:

- Concentrated or oxidizing acids
- Non-polar solvents, such as hydrocarbons and halogenated hydrocarbons
- Hydrofluoric acid (HF)
- Reagents that attack PP, PE, EPDM, and silicone
- Temperature of instrument and reagents must be between 15° and 40° C
- Piston/cylinder cartridge and pump housing are not autoclavable

### Items supplied with QuikSip™

- QuikSip™ controller with 200 cm of silicone tubing
- Two PP thread adapters 45/33 mm and 45/S40 mm
- Spare piston/cylinder cartridge
- Operating manual
- One-year warranty

See page 47 for other adapters

### Ordering information

Description	Cat. No.	2020 List Price
QuikSip™ bottletop aspirator	<a href="#">4723150</a>	\$341.20
8-channel manifold; PP, autoclavable, pack of 1	<a href="#">704526</a>	20.80
25 mL replacement piston/cylinder cartridge, pack of 3	<a href="#">704504</a>	90.40

### Bottles & adapters

See detailed listing on page 47.

