

# BRANDplates® Microplates

## Samples and plate selection

Given the broad range of available plates, view the Microplates at a Glance Chart on page 89 to select the optimal plate for your application. Samples are available without cost or obligation.

## Barcodes

BRAND® can provide high-quality and customized barcode labeling for your plate needs. Contact BrandTech® for more information.

## Accessories

BRAND® offers many accessories for use with microplate products such as lids, sealing films, and more. Lids with and without condensation rings are optimally matched to the corresponding BRANDplates® microplate to protect samples. For enhanced protection against contamination and evaporation, self-adhesive sealing films are available in a variety of materials and can be easily applied with a sealing paddle.

Condensation Rings	Height	Sterile	Pack of	Cat. No.	2020 List Price
<b>Lids for 96-well standard plates</b>					
Yes	8 mm	No	Pack of 100 (5 lids per bag)	<a href="#">782150</a>	\$170.00
No	8 mm	No	Pack of 100 (5 lids per bag)	<a href="#">782151</a>	172.00

For BRANDplates® microplates Cat. No. ranges: 781600-08, 781660-68, 781720-29, 781780-82, 781840-42, 781900-02, 781960-68, 782022-28, 782082



<b>Lids for 96-well plates with transparent bottom and all 384-well plates</b>					
No	4.5 mm	No	Pack of 50 lids (10 lids per bag)	<a href="#">782152</a>	112.40

For BRANDplates® microplates Cat. No. ranges: 781610-27, 781670-87, 781731-42, 781910-11, 781970-87, 782030-31, 782090-91

<b>Lids for all 1536-well plates</b>					
No	5.5 mm	No	Pack of 50 lids (10 lids per bag)	<a href="#">782153</a>	112.40

For BRANDplates® microplates Cat. No. ranges: 781640-42, 781700-02, 782000-02

Application	Description	Material	Pack of	Cat. No.	2020 List Price
<b>Self-adhesive Sealing Films</b>					
Real-Time PCR and crystallization	Single film, transparent, non-tacky	Polyolefin	100 sheets	<a href="#">781391ES</a>	\$260.40
PCR and ELISA	Single film	Polyester	100 sheets	<a href="#">781390</a>	131.20
<b>NEW!</b> PCR Sealing Film Strips, 1x8 Rows/Strip	8 strips/sheet	Polyester	400 sheets	<a href="#">781383</a>	45.80
ELISA	Single film	Polyester	100 sheets	<a href="#">701390ES</a>	78.80
ELISA	Mini-Strips (x4)	Polyester	50 4-strips	<a href="#">701391ES</a>	63.00
ELISA	Single film with assorted color tabs	Polyester	5 packs of 10 sheets	<a href="#">701392ES</a>	63.00
Cell and Tissue Culture	Single film, gas-permeable	Rayon	100 sheets	<a href="#">701364</a>	131.20
Cell and Tissue Culture	Single film, gas-permeable, sterile	Rayon	50 sheets	<a href="#">701365</a>	99.80
Pierceable and Resealable for Automation	Single film	Vinyl	100 sheets	<a href="#">701396ES</a>	136.40
Cold Storage	Single film	Aluminum	100 sheets	<a href="#">781381</a>	224.40
Cold Storage and Automation, Universal Fit	Single film	Aluminum	100 sheets	<a href="#">701395ES</a>	120.80
<b>NEW!</b> Aluminum Foil Strips, 1 x Rows/Strip	6 strips/sheet	Aluminum	300 sheets	<a href="#">781382</a>	56.00
DMSO Sample Storage	Single film, DMSO-resistant	Polypropylene	100 sheets	<a href="#">701367</a>	331.80
Adhesive-free Zones for Automation	Single film	PE/PP	50 sheets	<a href="#">701370</a>	157.40
Fluorescence Measurements	Single film	Black Vinyl	50 sheets	<a href="#">701371</a>	89.40
Luminescence Measurements	Single film	White Vinyl	50 sheets	<a href="#">701372</a>	73.40
Sealing Paddles for Sealing Films			Pack of 2	<a href="#">701381</a>	19.60



Sealing paddles