



VITsafe™ Safety Wash Bottles, LDPE, spray tube PP

- Durable safety imprint
- Material name in German, English, French, and Spanish
- Chemical formula, CAS No. hazard pictogram, signal word, risk phrases (H phrases), safety phrases (P phrases), as well as NFPA code

| Reagent | Cap Color | Volume mL | Diameter mm | Height* mm | Qty./Pack | Cat. No. | 2021 List Price |
|-----------------|-----------|-----------|-------------|------------|-----------|--------------------------|-----------------|
| Distilled Water | Blue | 500 | 76 | 166 | 12 | V1352819 | \$212.20 |
| Ethanol | Orange | 500 | 76 | 166 | 12 | V1352869 | 212.20 |
| Acetone | Red | 500 | 76 | 166 | 12 | V1452829 | 212.20 |
| Heptane | Red | 500 | 76 | 166 | 6 | V1352899 | 108.00 |
| Hexane | Red | 500 | 76 | 166 | 6 | V1352909 | 108.00 |
| Xylene | Red | 500 | 76 | 166 | 6 | V1352959 | 108.00 |

*Height without spray tube



Dropping Bottles, LDPE, with dropper screw cap, HDPE



| Volume mL | Thread GL | Diameter mm | Height mm | Qty./Pack | Cat. No. | 2021 List Price |
|-----------|-----------|-------------|-----------|-----------|------------------------|-----------------|
| 50 | 18 | 37 | 129 | 24 | V94587 | \$69.20 |
| 100 | 18 | 43 | 155 | 24 | V94687 | 77.80 |
| 250 | 25 | 58 | 183 | 12 | V94987 | 57.40 |
| 500 | 25 | 74 | 228 | 12 | V95087 | 75.60 |
| 1000 | 32 | 92 | 269 | 12 | V95187 | 102.40 |



Reagent Bottles, PP, with NS stoppers, PP

- Good chemical resistance, ideal for long-term storage of liquids
- Autoclavable at 121 °C according to DIN EN 285

| Volume mL | Stopper* NS | Height mm | Diameter mm | Qty./Pack | Cat. No. | 2021 List Price |
|-----------|-------------|-----------|-------------|-----------|-------------------------|-----------------|
| 100 | 14/23 | 106 | 52 | 20 | V100394 | \$73.40 |
| 1000 | 29/32 | 213 | 108 | 10 | V100694 | 98.00 |

*Stopper type A: With square knob and red core.



Reagent Bottles, PP, with screw caps, PP



- Screw caps help protect against outside contamination and leakage
- Autoclavable to 121 °C according to DIN EN 285, except for V100989 and V100889

| Volume mL | Thread GL | Height mm | Diameter mm | Qty./Pack | Cat. No. | 2021 List Price |
|-----------|-----------|-----------|-------------|-----------|--------------------------|-----------------|
| 100 | 18 | 100 | 52 | 20 | V100389 | \$65.00 |
| 250 | 25 | 132 | 70 | 20 | V100489 | 86.20 |
| 500 | 25 | 165 | 87 | 10 | V100589 | 62.00 |
| 1000 | 32 | 202 | 108 | 10 | V100689 | 82.00 |
| 2000 | 32 | 245 | 131 | 6 | V100789 | 102.40 |
| 5000 | 45 | 315 | 178 | 1 | V100889* | 29.80 |
| 10000 | 63 | 394 | 222 | 1 | V100989* | 74.60 |

*Bottle V100889 includes 1 handle; bottle V100989 includes 2 handles.

VITLAB® General Laboratory Products

Wide Mouth Reagent Bottles, PP, with screw caps, PP

- Screw caps help protect against outside contamination and leakage
- Wider mouth for easier filling
- Autoclavable at 121 °C according to DIN EN 285



| Volume mL | Thread GL | Height mm | Diameter mm | Qty./Pack | Cat. No. | 2021 List Price |
|-----------|-----------|-----------|-------------|-----------|-------------------------|-----------------|
| 100 | 32 | 96 | 55 | 20 | V101589 | \$86.20 |
| 250 | 45 | 132 | 73 | 20 | V101689 | 109.60 |
| 500 | 45 | 172 | 87 | 10 | V101789 | 70.20 |
| 1000 | 63 | 204 | 108 | 10 | V101889 | 126.80 |
| 2000 | 63 | 243 | 131 | 6 | V101989 | 95.80 |



Bottles, LDPE, with screw caps, PP

- Screw caps help protect against outside contamination and leakage



| | | | | | | |
|-----|----|-----|----|----|------------------------|-------|
| 50 | 18 | 85 | 37 | 24 | V94589 | 66.00 |
| 100 | 18 | 144 | 43 | 24 | V94689 | 70.20 |
| 250 | 25 | 135 | 65 | 12 | V94989 | 49.20 |
| 500 | 25 | 180 | 74 | 12 | V95089 | 70.20 |



Wide Mouth Bottles LDPE, with screw caps, PP

- Screw caps help protect against outside contamination and leakage
- Wider mouth for easier filling



| | | | | | | |
|------|----|-----|----|----|------------------------|--------|
| 250 | 45 | 146 | 58 | 12 | V93789 | 56.40 |
| 500 | 45 | 166 | 76 | 12 | V93989 | 77.80 |
| 1000 | 63 | 226 | 91 | 12 | V94189 | 110.60 |



Square Bottles, HDPE, with screw caps, LDPE

- Screw caps help protect against outside contamination and leakage
- Space-saving smaller footprint with square-shaped bottle

| Volume mL | Thread GL | Height mm | Size mm | Qty./Pack | Cat. No. | 2021 List Price |
|-----------|-----------|-----------|---------|-----------|------------------------|-----------------|
| 100 | 25 | 76 | 43 x 43 | 24 | V91789 | \$78.80 |
| 250 | 28 | 80 | 80 x 80 | 24 | V91989 | 105.60 |
| 500 | 32 | 106 | 90 x 90 | 12 | V92089 | 70.20 |
| 1000 | 32 | 187 | 80 x 80 | 12 | V92189 | 95.80 |



Wide Mouth Square Bottles, HDPE, with screw caps, LDPE

- Screw caps help protect against outside contamination and leakage
- Space-saving smaller footprint with square-shaped bottle

| | | | | | | |
|------|----|-----|---------|----|------------------------|--------|
| 100 | 32 | 78 | 46 x 46 | 24 | V92489 | 104.40 |
| 250 | 50 | 83 | 80 x 80 | 24 | V92689 | 118.20 |
| 500 | 65 | 120 | 90 x 90 | 12 | V92789 | 88.40 |
| 1000 | 65 | 168 | 90 x 90 | 12 | V92889 | 110.60 |



VITLAB® General Laboratory Products



Forceps

| Length mm | Qty./Pack | Cat. No. | 2021 List Price |
|-----------|-----------|----------|-----------------|
|-----------|-----------|----------|-----------------|

Forceps, PMP with pointed ends

| | | | |
|-----|----|------------------------|---------|
| 115 | 10 | V67895 | \$28.80 |
| 145 | 10 | V67995 | 34.20 |



Forceps, POM with rounded ends

| | | | |
|-----|---|------------------------|-------|
| 115 | 5 | V68099 | 13.40 |
| 145 | 5 | V68199 | 16.20 |
| 180 | 5 | V68299 | 30.40 |
| 250 | 5 | V68399 | 48.20 |



Scoops

| Volume mL | Length mm | Qty./Pack | Cat. No. | 2021 List Price |
|-----------|-----------|-----------|----------|-----------------|
|-----------|-----------|-----------|----------|-----------------|

Measuring Scoops, PP, volume marked on handle


| | | | | |
|------|-----|----|------------------------|--------|
| 2 | 60 | 12 | V39194 | \$6.60 |
| 5 | 82 | 12 | V39294 | 8.20 |
| 10 | 100 | 12 | V39394 | 10.60 |
| 25 | 135 | 12 | V39494 | 11.20 |
| 50 | 160 | 12 | V39594 | 21.40 |
| 100* | 200 | 12 | V39694 | 27.80 |
| 250 | 260 | 6 | V39794 | 35.20 |
| 500 | 315 | 6 | V39894 | 49.20 |
| 1000 | 385 | 6 | V39994 | 60.60 |

*Colors available – inquire



Measuring Scoops, HDPE

| | | | | |
|------|-----|----|------------------------|-------|
| 15 | 115 | 12 | V40093 | 11.80 |
| 25 | 135 | 12 | V40193 | 17.20 |
| 65 | 185 | 12 | V40293 | 22.40 |
| 110 | 215 | 12 | V40393 | 31.80 |
| 150 | 250 | 12 | V40493 | 37.40 |
| 350 | 310 | 6 | V40593 | 36.20 |
| 750 | 350 | 6 | V40693 | 59.60 |
| 1250 | 400 | 6 | V40793 | 83.20 |



Slide Staining Kit, POM

| Description | Size mm | Height mm | Qty./Pack | Cat. No. | 2021 List Price |
|-------------------------------|----------|-----------|-----------|------------------------|-----------------|
| Universal Slide Staining Kit* | 100 x 87 | 51 | 5 | 471800 | \$196.80 |
| Staining jar | 100 x 87 | 51 | 5 | 471500 | 161.20 |
| Staining rack | 91 x 79 | 38 | 5 | 471400 | 67.60 |

*Universal Kit 471800 includes one staining jar and a staining rack for 25 slides, 76 x 26 mm.

VITLAB 

Please see page 130 for description of various plastics.

VITLAB® General Laboratory Products

Funnels, PP

- Handle with loop for hanging

| Diameter mm | Approx vol. mL | Length mm | Stem Diameter mm | Qty./pack | Cat. No. | 2021 List Price |
|-------------|----------------|-----------|------------------|-----------|------------------------|-----------------|
| 30 | 5 | 45 | 1.5 | 24 | V40894 | \$21.40 |
| 30 | 6 | 47 | 4 | 24 | V41094 | 21.40 |
| 40 | 14 | 65 | 4 | 24 | V41194 | 25.60 |
| 50 | 32 | 85 | 7 | 24 | V41294 | 31.80 |
| 75 | 88 | 108 | 7 | 12 | V41394 | 20.20 |
| 100 | 222 | 155 | 8 | 12 | V41494 | 38.40 |
| 120 | 342 | 180 | 11 | 12 | V41594 | 48.20 |
| 150 | 817 | 220 | 15 | 12 | V41694 | 64.00 |



Funnels

Large Funnels, PP or HDPE*

- Integrated handle (except V41994)

| | | | | | | |
|-------------|-------|-----|----|---|------------------------|--------|
| 200 (PP) | 1300 | 200 | 22 | 6 | V41794 | 73.40 |
| 250 (PP) | 3200 | 260 | 30 | 6 | V41894 | 127.80 |
| 350 (PP) | 12500 | 440 | 35 | 1 | V41994 | 44.60 |
| 400 (HDPE)* | 12500 | 365 | 42 | 1 | V42294 | 114.00 |
| 430 (HDPE)* | 17500 | 420 | 37 | 1 | V42393 | 151.20 |

*HDPE Large Funnels V42294 and V42393 are not food-safe products.



Large funnels

Powder Funnels, PP

| | | | | | |
|-----|-----|----|----|------------------------|-------|
| 65 | 70 | 15 | 10 | V70794 | 27.40 |
| 80 | 75 | 21 | 10 | V70894 | 39.40 |
| 100 | 92 | 24 | 10 | V70994 | 41.60 |
| 120 | 105 | 27 | 10 | V71094 | 51.20 |
| 150 | 138 | 28 | 5 | V71194 | 33.00 |



Powder funnels

Sample Vials, PP, with snap caps, LDPE

| Volume mL | Height mm | Diameter mm | Qty./pack | Cat. No. | 2021 List Price |
|-----------|-----------|-------------|-----------|------------------------|-----------------|
| 5 | 25 | 20 | 25 | V68594 | \$49.20 |
| 18 | 57 | 22 | 25 | V68894 | 63.00 |
| 50 | 97 | 30 | 10 | V69194 | 32.60 |
| 160 | 110 | 50 | 10 | V69294 | 38.20 |



Sample Vials, PP, with screw caps, PP

| Volume mL | Height mm | Diameter mm | Thread GL | Qty./pack | Cat. No. | 2021 List Price |
|-----------|-----------|-------------|-----------|-----------|-------------------------|-----------------|
| 60 | 90 | 38 | 40 | 10 | V130394 | \$113.00 |
| 180 | 112 | 54 | 56 | 10 | V130594 | 165.20 |



Sample Vials, with snap caps

Sample Vials with screw caps

Please see page 130 for description of various plastics.



Bucket, PP



Bucket, HDPE

Buckets, PP, with handles, pouring lips, and graduations.



- With handle and lip for ease of pouring
- Markings at 1L graduations

| Volume L | Height mm | Diameter mm | Qty./Pack | Cat. No. | 2021 List Price |
|----------|-----------|-------------|-----------|------------------------|-----------------|
| 12 | 330 | 310 | 1 | V96694 | \$50.20 |
| 15 | 370 | 310 | 1 | V96794 | 56.40 |

Buckets, HDPE, with handles



- With carrying handle
- Markings at 1L graduations

| | | | | | |
|---------------------------|-----|-----|---|------------------------|-------|
| 5 | 240 | 250 | 1 | V96093 | 24.80 |
| 10 | 300 | 290 | 1 | V96393 | 29.80 |
| Lid for 5L bucket (LDPE) | — | — | 1 | V96293 | 6.60 |
| Lid for 10L bucket (LDPE) | — | — | 1 | V96593 | 7.80 |



Imhoff Sedimentation Cones, SAN

- Transparent
- Graduated to 1000 mL
- Screw closure allows drainage

| | | | |
|--|---|------------------------|--------|
| 1000 mL Cone | 1 | 388000 | 37.00 |
| Rack, PMMA, (for 2 cones) 150x300x290 mm | 1 | 388050 | 317.60 |

VITLAB® custom imprinting

VITLAB® has specialized capabilities in high-precision printing on plastic products with superior chemical and break resistance. Many VITLAB® products can be imprinted with logos or custom graduations for OEM or promotional use in relatively small quantities. Custom imprinted VITLAB® products are ideal for handling granulates, powders or liquids commonly used in agriculture, laboratories, and in the food, beverage and chemical industries. From multi-color logos to specialized measuring scales, VITLAB® can imprint plastic labware to meet your individual requirements and specifications. Custom imprinted labware by VITLAB® is of the highest quality, and offers a unique and functional means of presenting your company brand. Beakers, pitchers and scoops can also be made in colored plastic.

Contact BrandTech® Scientific for more information.

VITLAB 

