

# QuikSip™ Bottletop Aspirator

The QuikSip™ is a convenient bottletop aspirator that removes supernatants, cell culture media, and other biological liquids from dishes, flasks, culture bottles, centrifugation tubes, multiwell plates, and more. It is small, portable, and does not require an external vacuum or power source.

- **Moves easily between labs:** The QuikSip™ works without any external vacuum source, allowing it to be easily transferred between workstations. Each plunger stroke powers up to 25mL of aspiration.
- **Provides fingertip aspiration control:** The QuikSip™ includes a unique thumb-operated controller that is compatible with disposable pipette tips, micropipettes, and glass Pasteur pipettes.
- **Adapts to sterile environments:** All wetted parts are autoclavable (except piston/cylinder cartridge). Sterile piston/cylinders are available for use in sterile applications.
- **Fits most reagent containers:** Mounts on most common reagent bottles with its 45mm thread and included 45/33mm and 45/S40mm adapters. Other thread adapters available.

The economy aspirator for biological fluids

## Technical Data

### Operating limits and exclusions

Liquids unsuitable for use with the QuikSip™ bottletop aspirator include, but are not limited to:

- Concentrated or oxidizing acids.
- Non-polar solvents, such as hydrocarbons and halogenated hydrocarbons.
- Hydrofluoric acid (HF).
- Those that attack PP, PE, EPDM and silicone.
- Temperature of instrument and reagents must be between 15° and 40°C (59° and 104°F).
- Piston/cylinder cartridge and pump housing are not autoclavable.

### Items supplied with QuikSip™

- QuikSip™ controller with 200cm (78in.) of silicone tubing
- Two PP thread adapters 45/33mm and 45/S40mm. See page 41 for other adapters
- Spare piston/cylinder cartridge
- Operating manual
- One-year warranty



QuikSip™ aspirator shown with optional autoclavable 8-channel manifold.

### Ordering information

Description	Cat. No.	2016 List Price
QuikSip™ bottletop aspirator	4723150	\$325.00
8-channel manifold; PP, autoclavable, pack of 1	704526	19.11
25mL replacement piston/cylinder cartridge, pack of 3	704504	82.79

### Bottles & adapters

see detailed listing on page 41.