

HandyStep®

Repeating Pipettes



BRAND® HandyStep® touch/touch S

Choose the HandyStep® touch for your repetitive pipetting tasks—dispensing, auto-dispensing and the ability to store favorites. The HandyStep® touch S adds such advanced features like multi-aspiration, sequential dispensing and titration. Paired with the new PD-Tip™ // Dispenser Tips with tip size coding, they make the perfect team for easy pipetting.

- **Works with All Standard Tips:** The HandyStep® touch is compatible with BRAND® PD-Tip™ // Precision Dispenser Tips, Eppendorf® Combitips advanced®, and more.
- **Recognizes Tips Automatically:** The pipette automatically recognizes the size of encoded BRAND® PD-Tip™ // Precision Dispenser Tips upon tip insertion and displays the tip's volume, default dispensing volume, and default number of dispensing steps.
- **Reduces Repetitive Motion by up to 97%:** The instrument features a unique Auto-Dispense operating mode in which the pipette "learns" the operator's dispensing rhythm and dispenses at that rhythm as long as the control button is depressed – no programming needed!

Use the tips you want

The only touch screen repeating pipette compatible with all standard syringe tips, the HandyStep® touch accepts:

- BRAND® PD-Tip™ // Precision Dispenser Tips
- Rainin® Encode™ syringe tips
- Eppendorf® Combitips advanced®
- Fisherbrand® Dispenser Tips
- VWR® Combi-Syringes
- Gilson® Repet-Tips

Description	Cat. No.	2021 List Price
HandyStep® touch	705200	1,055.00
HandyStep® touch S	705210	1,299.00

BRAND® HandyStep® S

The HandyStep® S repeating pipette provides unprecedented versatility in a purely mechanical repeating pipette. In conjunction with BRAND® PD-Tip™ // Precision Dispenser Tips, the HandyStep® S offers nearly sixty different volume settings with accuracy, precision and dependability.

- **Comfortable:** The slim design of the HandyStep® S is lightweight, well-balanced, and provides a convenient control layout for smooth dispensing and one-handed volume selection.
- **Reliable:** Robust construction, inspired design and innovative polymers deliver a dependable instrument.
- **Ready-To-Go:** A purely mechanical design with no batteries to charge or replace means the HandyStep® S is ready when you are. A convenient integrated volume table provides a quick reference for dispensed volumes.
- **Flexible:** Ten sizes of BRAND® PD-Tip™ // Precision Dispenser Tips allow the HandyStep® S to dispense 59 different volumes. Open design allows use of third party tips—no more being locked into a sole source!
- **Contact-Free:** A tip ejection system allows the user to remove tips without handling.

Description	Cat. No.	2021 List Price
HandyStep® S	705110	\$402.80

BRANDTECH®
SCIENTIFIC, INC.

Toll Free: 888-522-2726
www.brandtech.com

PD-Tip™ // Syringe Tips

For Repeating Pipettes

The Perfect Tips for Your Repeating Pipette.

BRAND® encoded PD-Tip™ // syringe tips are manufactured from high-quality virgin plastics for use with all standard repeating pipettes in a wide variety of applications.

- **Compatible with Most Pipettes:** BRAND® PD-Tip™ // Precision Dispenser Tips are compatible with repeating pipettes using standard tips. BRAND® PD-Tip™ // Precision Dispenser Tips are encoded for automatic recognition by the BRAND® HandyStep® touch repeating pipette, the Rainin® AutoRep™ E pipette, and the Gilson® Repetman®.
- **Reliable:** A precision seal between the piston and cylinder ensure smooth operation, increased accuracy, and increased precision. Extended piston design minimizes dead volume and maximizes fluid utilization.
- **Manufactured from High Quality Materials:** PD-Tip™ // syringe tips are made entirely from high-quality virgin polyethylene and polypropylene. (Piston of 0.1 mL tip is liquid crystal polymer-LCP).
- **BIO-CERT® PD-Tip™ // Precision Dispenser Tips** are sterile, individually wrapped, and certified free of DNA, RNase, endotoxins, and ATP to meet the most stringent sample integrity requirements.
- **Batch Certificate of Performance** which documents tip accuracy and precision testing is supplied with each pack of BRAND® PD-Tip™ // Precision Dispenser Tips.



PD-Tip™ // Precision Dispenser Tips BIO-CERT® sterile, individually wrapped

Capacity,	Pack of	Cat. No.	2021 List Price
0.1 mL	100	705730	\$171.60
0.5 mL	100	705732	171.60
1 mL	100	705734	171.60
1.25 mL	100	705736	171.60
2.5 mL	100	705738	171.60
5 mL	100	705740	171.60
10 mL	100	705742	171.60
12.5 mL	100	705744	171.60
25 mL (one adapter included)	25	705746	52.60
50 mL (one adapter included)	25	705748	52.60
Adapter for 25 mL and 50 mL , BIO-CERT® sterile, PP, autoclavable			
25 mL and 50 mL Adapter	5	702399	78.40

*The Eppendorf® Repeater® Plus, Repeater® M4, Repeater® Pro, Stream, and XStream are designed to preclude their operation with any tips other than Eppendorf® Combitips® advanced.

PD-Tip™ // Precision Dispenser Tips, non-sterile

Capacity,	Pack of	Cat. No.	2021 List Price
0.1 mL	100	705700	\$125.20
0.5 mL	100	705702	125.20
1 mL	100	705704	125.20
1.25 mL	100	705706	125.20
2.5 mL	100	705708	125.20
5 mL	100	705710	125.20
10 mL	100	705712	125.20
12.5 mL	100	705714	125.20
25 mL (one adapter included)	50	705716	84.80
50 mL (one adapter included)	25	705718	42.40
Multi-pack (20 each 0.5, 1, 1.25, 2.5, 5, 10 and 12.5 mL)			
Adapter for 25 mL and 50 mL , non-sterile, PP, autoclavable	140	705720	169.60
25 mL and 50 mL	10	702398	99.60



Lab Rats Trust BrandTech!

Toll Free: 888-522-2726 www.brandtech.com

HandyStep_BT_Fly-2021

BRANDTECH®
SCIENTIFIC, INC.