



# BRAND<sup>®</sup> PCR PLASTICS

## Comprehensive Line of High Quality PCR Plastics

BRAND<sup>®</sup> offers a collection of PCR plastics produced from pure raw material without lubricants and chemicals using the most advanced cleanroom techniques. In addition to the standard line of tubes, strips and plates, BRAND<sup>®</sup> is offering high quality plates with advanced features for fast, reliable and reproducible PCR results.


- Standard (0.2mL) and low profile (0.15mL) plates and tubes compatible with most thermal cyclers
- PCR products are free from DNA, DNase, RNase and other PCR inhibitors
- Thin uniform tube wall thickness and rigid upper plate for optimal thermal transfer, short cycle times, stable product handling and consistency
- Unique blue alphanumeric coding and cut corner marking for easy identification and fast orientation
- Innovative sealing films and TPE sealing mats drastically reduce sample evaporation
- Cap tool for reliable, soft-touch sealing and detaching of PCR caps
- BRAND<sup>®</sup>'s PCR box/rack securely holds 96-well PCR plates, tube strips, and individual tubes for easy handling, and withstands temperatures from -80° C to +121° C
- PCR Mini Cooler helps protect samples from warming. Insulating gel maintains 4° C temperature for approximately 3 hours and changes color from violet to pink when temperature reaches 7° C




BRAND PCR tube strips and caps



BRAND PCR tubes



BRAND PCR plates with blue alphanumeric coding



Low profile plates available



PCR Mini-Cooler



PCR Box/Rack




Cap Tool

**White PCR** products are designed specifically to optimize results when performing quantitative Real Time PCR (qPCR). White PCR strips and plates offer significantly better results in the analysis of fluorescence signals and consequently with the quantification of the replicated DNA.

- Uniformly colored with titanium dioxide (TiO<sub>2</sub>)
- Smooth surfaces to optimize reflection of the fluorescence signal
- Low profile plates and strips



White PCR strips optimized for qPCR



White 96-well PCR plates optimized for qPCR

Ordering information on reverse.



**BRANDTECH<sup>®</sup>**  
SCIENTIFIC, INC.

# BRAND<sup>®</sup> PCR Tubes, Strips and Plates

Description	Qty./ Pack	Cat. No.	2018 List Price
<b>PCR Tubes</b>			
0.2mL tubes with domed cap, Clear	1000	781300	\$107.00
0.2mL tubes with flat cap, Clear	1000	781305	57.00
0.2mL tubes with flat cap, Rose	1000	781301	57.00
0.2mL tubes with flat cap, Yellow	1000	781302	57.00
0.2mL tubes with flat cap, Green	1000	781303	57.00
0.2mL tubes with flat cap, Blue	1000	781304	57.00
0.5mL tubes with flat cap, Clear	1000	781310	60.00
0.5mL tubes with flat cap, Rose	1000	781311	60.00
0.5mL tubes with flat cap, Yellow	1000	781312	60.00
0.5mL tubes with flat cap, Green	1000	781313	60.00
0.5mL tubes with flat cap, Blue	1000	781314	60.00
<b>PCR 8-Strip Tubes</b>			
8 x 0.2mL flexible tube strip with individual attached flat caps, clear	120	781332	95.00
8 x 0.2mL flexible tube strip with individual attached domed caps, clear <b>NEW!</b>	120	P1300	90.00
8 x 0.15mL flexible tube strip with ind. attached flat caps, clear	120	781333	95.00
8 x 0.2mL rigid tube strip with ind. attached flat caps, clear <b>NEW!</b>	120	P1200	90.00
8 x 0.15mL rigid tube strip with ind. attached flat caps, clear <b>NEW!</b>	120	P1100	90.00
8 x 0.2mL tube strip with 8-strip flat caps, Clear (qPCR)	250	781326	251.00
8 x 0.2mL tube strip with 8-strip domed caps, Clear	250	781327	207.00
8 x 0.2mL tube strip, Clear (qPCR)	125	781320	71.00
8 x 0.2mL tube strip, Rose	125	781321	71.00
8 x 0.2mL tube strip, Yellow	125	781322	71.00
8 x 0.2mL tube strip, Green	125	781323	71.00
8 x 0.2mL tube strip, Blue	125	781324	71.00
<b>PCR 8-Strip Caps</b>			
8-strip flat caps, Clear (qPCR)	125	781334	58.00
8-strip domed caps, Clear	125	781340	34.80
8-strip domed caps, Rose	125	781341	34.80
8-strip domed caps, Yellow	125	781342	34.80
8-strip domed caps, Green	125	781343	34.80
8-strip domed caps, Blue	125	781344	34.80
8-strip flat caps, Clear (qPCR), for plates <b>NEW!</b>	300	781413	138.00
8-strip domed caps, Clear, for plates <b>NEW!</b>	300	781414	138.00
<b>PCR 12-Strip Tubes</b>			
12 x 0.2mL tube strip, Clear (qPCR)	125	781280	206.00
12 x 0.2mL tube strip, Rose	125	781281	206.00
12 x 0.2mL tube strip, Yellow	125	781282	206.00
12 x 0.2mL tube strip, Green	125	781283	206.00
12 x 0.2mL tube strip, Blue	125	781284	206.00
<b>PCR 12-Strip Caps</b>			
12-strip domed caps, Clear	125	781290	53.00
12-strip domed caps, Rose	125	781291	53.00
12-strip domed caps, Yellow	125	781292	53.00
12-strip domed caps, Green	125	781293	53.00
12-strip domed caps, Blue	125	781294	53.00

Description	Qty./ Pack	Cat. No.	2018 List Price
<b>Clear PCR Plates</b>			
24-well plate, non-skirted, 0.2mL	40	781411	\$96.00
48-well plate, non-skirted, 0.2mL	20	781415	96.00
96-well plate, non-skirted, elevated rim, 0.2mL	50	781350	195.00
96-well plate, non-skirted, low-profile, 0.15mL	50	781366	196.15
96-well plate, non-skirted, standard profile, 0.2mL	50	781368	196.15
96-well plate, semi-skirted, elevated rim, 0.2mL	50	781400	241.00
96-well plate, semi-skirted, low-profile, 0.15mL	50	781371	196.15
96-well plate, semi-skirted, low-profile, raised skirt, 0.15mL	50	781373	196.15
96-well plate, semi-skirted, standard profile, 0.2mL	50	781375	196.15
96-well plate, skirted, low-profile, 0.15mL	50	781377	346.25
384-well plate, skirted	50	781345	294.00
384-well plate, full skirt	50	781347	371.00
384-well plate, full skirt, rigid plate	50	781348	538.00
<b>White PCR Plastics (all suitable for qPCR)</b>			
8 x 0.2mL tube strip	125	781325	158.00
8 x 0.2mL tube strip with individual attached flat caps	120	781316	95.00
8 x 0.15mL low-profile tube strip with individual attached flat caps	120	781317	95.00
24-well plate, non-skirted, 0.2mL	40	781412	110.00
48-well plate, non-skirted, 0.2mL	20	781416	110.00
96-well plate, non-skirted, 0.2mL	50	781354	242.00
96-well plate, non-skirted, low-profile, 0.15mL	50	781367	222.20
96-well plate, non-skirted, standard profile, 0.2mL	50	781369	222.20
96-well plate, semi-skirted, standard profile, 0.2mL	50	781376	222.20
96-well plate, semi-skirted, standard profile, 0.2mL	50	781357	328.00
96-well plate, semi-skirted, low-profile, 0.15mL	50	781372	222.20
96-well plate, semi-skirted, low-profile, raised skirt, 0.15mL	50	781374	222.20
96-well plate, skirted, low-profile, 0.15mL	50	781378	360.25
96-well plate, low-profile, for Roche LightCycler <sup>®</sup> , 0.15mL	50	781364	296.00
96-well plate, low-profile, for Roche LightCycler <sup>®</sup> , with fifty sheets of qPCR sealing film, 0.15mL	50	781365	390.00
384-well plate, full skirt, for Roche LightCycler <sup>®</sup>	50	781358	562.00
<b>Accessories</b>			
Cap Tool	1	781419	56.12
PCR Box/Rack	5	781362	77.00
PCR Mini Cooler	2	781260	128.00
<b>PCR Mats and Sealing Films</b>			
TPE 96-well mat for 781350 and 781400	5	781405	51.00
TPE 24-well mat for 781411	10	781402	29.00
TPE 48-well mat for 781415	10	781403	53.00
PP Film for PCR, single sheet	100	781390	125.00
Polyolefin film for Real-Time PCR (qPCR), single sheet	100	781391ES	248.00
Sealing Paddle for Sealing Films <b>NEW!</b>	2	701381	18.65



Toll Free 888-522-2726  
[www.brandtech.com](http://www.brandtech.com)