

Reliable Tools

for Craft Brewing Labs



Faster, Easier Titrations with Class A Precision

Titrette® bottle-top burette



- Meets Class A standards for glass burettes!
- Popular for sulfur dioxide (SiO₂) analyses, acid-base titrations, solvent analyses and more
- Compatible with NaOH, KOH, KI, H₂SO₄, KMnO₄, Iodine and more
- Eliminates meniscus reading errors
- Models available with direct data transmission to PC



Description	Cat. No.	2018 List Price
10mL (standard)	4760141	\$1,025.00
25mL (standard)	4760151	1,025.00
50mL (standard)	4760161	1,025.00

Comfortable, Convenient Meniscus Control

macro™ Pipette Controller

- Replaces old-fashioned pipette bulbs with simple, one-finger control
- Rapid or sensitive filling; fast or drop-wise dispensing
- Fits most pipettes 0.1-200mL

Description	Cat. No.	2018 List Price
macro™ pipette controller		
Gray	26200	\$84.00
Blue	26202	84.00
Magenta	26203	84.00
Green	26201	84.00



The Standard in Bottletop Dispensers

Dispensette® S dispensers

- Rugged: no seals to wear out!
- Compatible with most laboratory reagents
- Simple to clean & calibrate
- Recirculation valve reduces splashing and helps conserve reagents
- Autoclavable to 250°F (121°C)
- Adapters for common bottle threads included
- Optional flexible discharge tube simplifies serial dispensing



Dispensette® S Volume, mL	Without Recirculation Valve		With Recirculation Valve	
	Cat. No.	2018 List Price	Cat. No.	2018 List Price
Dispensette® S, Digital				
0.1-1	4600310	\$502.00	4600311	\$533.00
0.2-2	4600320	502.00	4600321	533.00
0.5-5	4600330	502.00	4600331	533.00
1-10	4600340	502.00	4600341	533.00
2.5-25	4600350	678.00	4600351	709.00
5-50	4600360	688.00	4600361	719.00
Dispensette® S, Analog-adjustable				
0.1-1	4600100	\$440.00	4600101	\$460.00
0.2-2	4600120	440.00	4600121	460.00
0.5-5	4600130	440.00	4600131	460.00
1-10	4600140	440.00	4600141	460.00
2.5-25	4600150	616.00	4600151	636.00
5-50	4600160	631.00	4600161	652.00
10-100	4600170	983.00	4600171	1,004.00

Fixed volume models available. For strong acids and nonpolar organics, ask for the BRAND® Dispensette® S Organic dispensers.

food_bev_Enology_bt_2018



Lab Rats Trust BrandTech!

BRANDTECH®
SCIENTIFIC, INC.

Toll Free: 888-522-2726

www.brandtech.com

Reliable Tools for **Craft Brewing Labs**

Precise, Slim, Light & Dependable Repeater

HandyStep® S Repeating Pipette



The HandyStep® S provides unprecedented versatility in a mechanical repeating pipette. It offers 59 different volume settings with accuracy, precision and dependability.

- Comfortable, slim design is light weight, well balanced and offers easy one handed volume selection
- Flexible; accepts third party tips – no more being locked into a sole source
- Tip ejection system allows tip removal without user handling

HandyStep® electronic also available

Description	2018 Cat. No.	2018 List Price
HandyStep® S Repeating Pipette	705110	\$378.00
HandyStep® electronic Repeating Pipette	705002	1,019.00
PD-Tip™ size	Non-sterile tips	BIO-CERT® Sterile
0.1mL (100 tips)	702402 \$125.00	702683 \$167.00
0.5mL (100 tips)	702370 117.00	702684 160.00
1.0 mL (100 tips)	702406 117.00	702685 160.00
1.25mL (100 tips)	702372 117.00	702686 160.00
2.5mL (100 tips)	702374 117.00	702688 160.00
5mL (100 tips)	702376 117.00	702690 160.00
10.0 mL (100 tips)	702407 117.00	702691 160.00
12.5mL (100 tips)	702378 117.00	702692 160.00
25mL (50 tips)	702380 73.00	—
25mL (25 tips)	—	702694 47.00
50mL (25 tips)	702382 42.00	702696 48.00
Multi-pack (20 tips each 0.5,1, 1.25, 2.5, 5, 10 & 12.5mL)	702368 158.00	—

*BIO-CERT® is certified sterile, and DNA-, RNase-, ATP- and endotoxin-free.

Comfortable, Convenient Meniscus Control

accu-jet® pro pipette controller

- Provides smooth, intuitive meniscus control
- Simple, sensitive soft blow-out
- Contoured grip for comfort
- Mode switch permits either powered or gravity dispensing

Description	2018 Cat. No.	2018 List Price
accu-jet® pro pipette controller		
Dark Blue	26330	\$407.00
Magenta	26331	407.00
Green	26332	407.00
Royal Blue	26333	407.00



Lightweight and Rugged

Transferpette® S Adjustable Volume Single Channel Pipettes



Capacity, µL	2018 Cat. No.	2018 List Price
0.1-1	704768	\$342.00
0.1-2.5	704769	311.00
0.5-10	704770	311.00
2-20	704772	311.00
5-50	704773	311.00
10-100	704774	311.00
20-200	704778	311.00
100-1000	704780	311.00
500-5000 (0.5-5mL)	704782	311.00
1000-10000 (1-10mL)	704784	311.00

- Simple, thumb-tip volume adjustment, even with gloves!
- Shortest stroke length among leading pipettes reduces strain
- Fully and repeatedly autoclavable
- Lightweight and rugged
- Accepts universal tips for economy
- Fixed volume single channels also available



Transferpette® S -8/-12 Multichannel Pipettes

Capacity, µL	2018 8 Channel Cat. No.	2018 8 Channel List Price	2018 12 Channel Cat. No.	2018 12 Channel List Price
0.5-10	703700	\$760.00	703720	\$938.00
5-50	703706	760.00	703726	938.00
10-100	703708	760.00	703728	938.00
20-200	703710	760.00	703730	938.00
30-300	703712	760.00	703732	938.00

Transferpette® S -8/-12 pipettes supplied with performance certificate, operating manual, shelf mount, two tip boxes, reagent reservoir, pack of replacement O-rings. One-year warranty.



Compact Vacuum Pump for Filtration

VACUUBRAND® ME 1 Vacuum Pump



VACUUBRAND® ME 1 is a compact, high performance, easy to use vacuum pump which is ideal for single and multiple filtrations. An optional vent control valve with dial gauge is available for variable fine adjustment of pumping speed.

- Small footprint conserves valuable bench space
- Whisper quiet & low vibration minimizes background noise in the lab
- PTFE diaphragms and valves offer good chemical resistance
- High performance 14lpm pumping speed & 100mbar ultimate vacuum

Model	2018 Cat. No.	2018 List Price
ME 1 Vacuum Pump	721003	\$868.00
Vacuum regulation valve with manometer for ME 1	696842	182.08

Need Lab Vacuum?

Use the [BrandTech® Vacuum Selection Guide](#) to find the best vacuum pump for a wide range of applications! Go to www.brandtech.com.