

# TIP SELECTION GUIDE

Tip volume is less than pipette's nominal volume.  
\*Calibration and possible adjustment needed.

		Pipette tips (no filter)			
		Standard		ULR	
		Non-sterile	Sterile	Non-sterile	Sterile
Pack type	Cat. No.	Cat. No.	Cat. No.	Cat. No.	
<b>For pipettes:</b>		<b>0.1-20 µl</b>			
705868, 705869*, 705870*	Bulk	<a href="#">732002</a>	-	-	-
705808*	Bulk, xxl	<a href="#">732022</a>	-	-	-
705900, 705920	Tip-Refill	<a href="#">732102</a>	<a href="#">732122</a>	-	-
705319*, 705320*, 705339*, 705340*	TipBox	<a href="#">732202</a>	-	<a href="#">732302</a>	-
705419, 705420, 705469, 705470	TipBox, Sterile	-	<a href="#">732222</a>	-	<a href="#">732322</a>
<b>For pipettes:</b>		<b>0.5-20 µl</b>			
705869, 705870	Bulk	<a href="#">732004</a>	-	-	-
705808	Bulk, xxl	<a href="#">732024</a>	-	-	-
705900, 705920	Tip-Refill	<a href="#">732104</a>	<a href="#">732124</a>	-	-
705319, 705320, 705339, 705340	TipBox	<a href="#">732204</a>	-	<a href="#">732304</a>	-
705419, 705420, 705469, 705470	TipBox, Sterile	-	<a href="#">732224</a>	-	<a href="#">732324</a>
	TipStack	<a href="#">732244</a>	<a href="#">732264</a>	<a href="#">732344</a>	<a href="#">732364</a>
<b>For pipettes:</b>		<b>1-50 µl</b>			
705869*, 705870	Bulk	<a href="#">732006</a>	-	-	-
705808	Bulk, xxl	<a href="#">732026</a>	-	-	-
705900, 705920	Tip-Refill	<a href="#">732106</a>	<a href="#">732126</a>	-	-
705319, 705320, 705339, 705340	TipBox	<a href="#">732206</a>	-	<a href="#">732306</a>	-
705419, 705420, 705469, 705470	TipBox, Sterile	-	<a href="#">732226</a>	-	<a href="#">732326</a>
<b>For pipettes:</b>		<b>2-200 µl</b>			
705872, 705873, 705874, 705878	Bulk	<a href="#">732008</a>	-	-	-
705816, 705820, 705828, 705838, 705844	Bulk, xxl	<a href="#">732028</a>	-	-	-
705906, 705926, 705908, 705928, 705910, 705930, 705912, 705932	Tip-Refill	<a href="#">732108</a>	<a href="#">732128</a>	-	-
705323, 705343	TipBox	<a href="#">732208</a>	-	<a href="#">732308</a>	-
705423, 705473, 705424, 705474, 705426, 705476	TipBox, Sterile	-	<a href="#">732228</a>	-	<a href="#">732328</a>
	TipStack	<a href="#">732248</a>	<a href="#">732268</a>	<a href="#">732348</a>	<a href="#">732368</a>
<b>For pipettes:</b>		<b>5-300 µl</b>			
705872, 705873, 705874, 705878	Bulk	<a href="#">732010</a>	-	-	-
705816, 705820, 705828, 705838, 705844	Bulk, xxl	<a href="#">732030</a>	-	-	-
705906, 705926, 705908, 705928, 705910, 705930, 705912, 705932	Tip-Refill	<a href="#">732110</a>	<a href="#">732130</a>	-	-
705323, 705343	TipBox	<a href="#">732210</a>	-	<a href="#">732310</a>	-
705423, 705473, 705424, 705474, 705426, 705476	TipBox, Sterile	-	<a href="#">732230</a>	-	<a href="#">732330</a>
<b>For pipettes:</b>		<b>50-1000 µl</b>			
705880	Bulk	<a href="#">732012</a>	-	-	-
705854, 705862	Bulk, xxl	<a href="#">732032</a>	-	-	-
705326, 705346	Tip-Refill	<a href="#">732112</a>	<a href="#">732132</a>	-	-
	TipBox	<a href="#">732212</a>	-	<a href="#">732312</a>	-
	TipBox, Sterile	-	<a href="#">732232</a>	-	<a href="#">732332</a>
	TipStack	<a href="#">732252</a>	<a href="#">732272</a>	<a href="#">732352</a>	<a href="#">732372</a>
<b>For pipettes:</b>		<b>0.5-5 mL</b>			
705882	Bulk	<a href="#">702595</a>	-	-	-
705327, 705347	Bulk, xxl	<a href="#">702600</a>	-	-	-
	TipBox, 5 mL	<a href="#">702605</a>	-	-	-
<b>For pipette:</b>		<b>1-10 mL</b>			
705884	Bulk	<a href="#">702603</a>	-	-	-
	Bulk, xxl	<a href="#">702604</a>	-	-	-
	TipBox, 10 mL	<a href="#">702608</a>	-	-	-

# TIP SELECTION GUIDE

Tip volume is less than pipette's nominal volume.  
\*Calibration and possible adjustment needed.

		Filter Pipette tips			
		Standard		ULR	
		Non-sterile	Sterile	Non-sterile	Sterile
Pack type		Cat. No.	Cat. No.	Cat. No.	Cat. No.
<b>For pipettes:</b>		<b>0.1-1 µl</b>			
705868, 705869*, 705870*	Tip-Refill	<a href="#">732602</a>	<a href="#">732622</a>	-	-
705808*	TipBox	<a href="#">732702</a>	-	<a href="#">732802</a>	-
705900, 705920	TipBox, Sterile	-	<a href="#">732722</a>	-	<a href="#">732822</a>
705319*, 705339*	-	-	-	-	-
705419, 705469	-	-	-	-	-
<b>For pipettes:</b>		<b>0.5-10 µl</b>			
705869, 705870	Tip-Refill	<a href="#">732604</a>	<a href="#">732624</a>	-	-
705808	TipBox	<a href="#">732704</a>	-	<a href="#">732804</a>	-
705900, 705920	TipBox, Sterile	-	<a href="#">732724</a>	-	<a href="#">732824</a>
705319, 705320, 705339, 705340	-	-	-	-	-
705419, 705420, 705469, 705470	-	-	-	-	-
<b>For pipettes:</b>		<b>1-20 µl</b>			
705870	Tip-Refill	<a href="#">732606</a>	<a href="#">732626</a>	-	-
705808	TipBox	<a href="#">732706</a>	-	<a href="#">732806</a>	-
705900, 705920	TipBox, Sterile	-	<a href="#">732726</a>	-	<a href="#">732826</a>
705319, 705320, 705339, 705340	-	-	-	-	-
705419, 705420, 705469, 705470	-	-	-	-	-
<b>For pipettes:</b>		<b>2-20 µl</b>			
705872, 705873, 705874, 705878	Tip-Refill	<a href="#">732608</a>	<a href="#">732628</a>	-	-
705816, 705820, 705828, 705838, 705844	TipBox	<a href="#">732708</a>	-	<a href="#">732808</a>	-
705906, 705926, 705908, 705928, 705910, 705930	TipBox, Sterile	-	<a href="#">732728</a>	-	<a href="#">732828</a>
705323, 705343	-	-	-	-	-
705423, 705473, 705424, 705474	-	-	-	-	-
<b>For pipettes:</b>		<b>5-100 µl</b>			
705872, 705873, 705874, 705878	Tip-Refill	<a href="#">732610</a>	<a href="#">732630</a>	-	-
705816, 705820, 705828, 705838, 705844	TipBox	<a href="#">732710</a>	-	<a href="#">732810</a>	-
705906, 705926, 705908, 705928, 705910, 705930, 705912, 705932, 705323, 705343	TipBox, Sterile	-	<a href="#">732730</a>	-	<a href="#">732830</a>
705423, 705473, 705424, 705474, 705426, 705476	-	-	-	-	-
<b>For pipettes:</b>		<b>5-200 µl</b>			
705878, 705844	Tip-Refill	<a href="#">732612</a>	<a href="#">732632</a>	-	-
705906, 705926, 705908, 705928, 705910, 705930, 705912, 705932	TipBox	<a href="#">732712</a>	-	<a href="#">732812</a>	-
705323, 705343	TipBox, Sterile	-	<a href="#">732732</a>	-	<a href="#">732832</a>
705423, 705473, 705424, 705474, 705426, 705476	-	-	-	-	-
<b>For pipettes:</b>		<b>50-1000 µl</b>			
705880	Tip-Refill	<a href="#">732614</a>	<a href="#">732634</a>	-	-
705854, 705862	TipBox	<a href="#">732714</a>	-	<a href="#">732814</a>	-
705326, 705346	TipBox, Sterile	-	<a href="#">732734</a>	-	<a href="#">732834</a>