



Reliable Tools for **High Throughput Screening Labs**

High throughput screening labs have many challenges to meet the demanding research requirements of drug discovery, biomedical and biochemical research initiatives, including increasing screening capacity while decreasing screening cost. With all of the high tech equipment you have in your HTS facility—why would you want anything less than the highest quality assay plates?

Full Range of 96-, 384- and 1536-well Plates

- **Surface Treatments** – Available in non-treated, sterile and non-sterile, or modified for immunoassay and cell culture applications
- **Colors, volumes, and well shapes** – Transparent, White, Black, Colored with clear bottoms, and all bottom shapes, round (U), flat (F), conical (V), and curved (C).

Broad Instrument Compatibility

- **Manufactured according to ANSI-/SLAS – Standards** to ensure compatibility for instrumentation including liquid handlers, plate washers, plate readers, and more
- Dimensional drawings available upon request

Quality Material

- Premium quality molding of high quality virgin resins for the highest product performance
- Internal quality controls and independent accredited laboratory testing to ensure plates are DNA, DNase, RNase, endotoxin, and cytotoxin free

Special Features

- **Color Coded Embossing** – For easy identification of surface type
- **Custom Bar Coding Available** – For plate identification with robotic equipment or for QC processes
- **Custom Packaging Available**
- **Microplate Selection Guide** – Find the best microplates for your application at www.brandtech.com/plateselection
- **Sample Program** – Give our plates a try, FREE samples!



Toll Free: 888-522-2726
www.brandtech.com

BRANDplates® Microplates from BrandTech Scientific



BRANDplates® microplates at a Glance

96-well Microplates

Type		Non-treated		Immunological surfaces			Cell culture surfaces			
Standard										
Color	Bottom /Well volume	pureGrade™	pureGrade™ S	immunoGrade™	hydroGrade™	lipoGrade™	cellGrade™	cellGrade™ plus	cellGrade™ premium	inertGrade™
transparent	U / 330µl	781600	781660	781720	781780	781840	781960	–	–	781900
transparent	V / 360µl	781601	781661	781721	–	781841	781961	–	–	–
transparent	F / 350µl	781602	781662	781722	781782	781842	781962	782022	782082	781902
transparent	C / 350µl	781603	781663	781723	781783	781843	781963	–	–	–
white	U / 330µl	781604	–	781724	781784	781844	–	–	–	781904
white	F / 350µl	781605	781665	–	–	–	781965	782025	782085	781905
white	C / 350µl	–	–	781726	781786	781846	–	–	–	781906
black	U / 330µl	781607	–	781727	781787	781847	–	–	–	781907
black	F / 350µl	781608	781668	–	–	–	781968	782028	782088	781908
black	C / 350µl	–	–	781729	781789	781849	–	–	–	781909
With transparent bottom										
Color	Bottom /Well volume	pureGrade™	pureGrade™ S	immunoGrade™	hydroGrade™	lipoGrade™	cellGrade™	cellGrade™ plus	cellGrade™ premium	inertGrade™
white	F / 330µl	781610	781670	–	–	781850	781970	782030	782090	781910
black	F / 330µl	781611	781671	781731	–	781851	781971	782031	782091	781911
Strip plates										
Color	Bottom /Well volume	pureGrade™	pureGrade™ S	immunoGrade™	hydroGrade™	lipoGrade™	cellGrade™	cellGrade™ plus	cellGrade™ premium	inertGrade™
transparent, w/o grid	F / 360µl	782300	–	782305	–	–	–	–	–	–
transparent, w/ grid	F / 350µl	782301	–	782306	–	–	–	–	–	–

384-well HTS Microplates

Type		Non-treated		Immunological surfaces			Cell culture surfaces			
Standard										
Color	Bottom /Well volume	pureGrade™	pureGrade™ S	immunoGrade™	hydroGrade™	lipoGrade™	cellGrade™	cellGrade™ plus	cellGrade™ premium	inertGrade™
transparent	F / 100µl	781620	781680	781740	781800	781860	781980	782040	782100	–
white	F / 100µl	781621	781681	781741	781801	781861	781981	–	–	–
black	F / 100µl	781622	781682	781742	781802	781862	781982	–	–	–
With transparent bottom										
Color	Bottom / Well volume	pureGrade™	pureGrade™ S	immunoGrade™	hydroGrade™	lipoGrade™	cellGrade™	cellGrade™ plus	cellGrade™ premium	inertGrade™
White	F / 120µl	781626	781686	–	–	781866	781986	782046	782106	–
Black	F / 120µl	781627	781687	–	–	781867	781987	782047	782107	–

1536-well UHTS Microplates

Type		Non-treated		Immunological surfaces			Cell culture surfaces			
Standard										
Color	Bottom /Well volume	pureGrade™	pureGrade™ S	immunoGrade™	hydroGrade™	lipoGrade™	cellGrade™	cellGrade™ plus	cellGrade™ premium	inertGrade™
transparent	F / 10µl	781640	781700	–	–	–	782000	–	–	–
white	F / 10µl	781641	781701	–	–	–	782001	–	–	–
black	F / 10µl	781642	781702	–	–	–	782002	–	–	–

Interactive Selection Guide with more information about BRANDplates® microplates available at www.brandtech.com.

BT_fly_plates_8-15



Lab Rats Trust BrandTech®!



Toll Free: 888-522-2726

www.brandtech.com