

The **NEW** Transferpette® S

A perfect fit for every hand!

BRAND. For lab. For life.®



Precision workmanship, high grade materials, "Made in Germany" quality – it's what you see on first glance and appreciate every time you use the new Transferpette® S. The finer points of the design will be evident when the instrument is in your hand. Experience how simple, easy and efficient pipetting can be with the new Transferpette® S. Single- or multi-channel, there's a Transferpette® S pipette for your application.

- **Slim hand grip**—Provides a comfortable and secure fit in any hand
- **One-handed operation**—Perform all functions (set volume, pipette, tip ejection) without changing your grip
- **Volume-change protection**—Simple to operate and clearly shows when the volume setting is locked
- **Slim shaft**—Narrow tubes? With the slim shaft, easily pipette into narrow tubes such as 15 mL conical tubes
- **Easy cleaning**—The shaft tucks inside the handgrip for smooth surfaces, making cleaning easier
- **Low forces**—Short pipetting stroke, easy mounting and tip ejection make pipetting a breeze
- **Universal tip cone** fits tips from most quality manufacturers
- **Completely autoclavable** at 121 °C
- **Color-coded pipette** accents make pipette and tip selection easy
- **Easy mounting and ejection** of tips with O-rings on multi-channel pipettes
- **3-year warranty** for single-channel pipettes

Reaktionsgefäßständer
Micro Test Tube Rack



BRANDTECH®
SCIENTIFIC, INC.

Transferpette® S



Transferpette® S pipettes, adjustable volume (fixed volume also available)

| | Capacity, μ l (color-coded) | Description | Cat. No. | 2021 Price |
|---|------------------------------------|-------------|----------|------------|
| ● | 0.1–1 | D-1 | 705868 | \$297.00 |
| ● | 0.1–2.5 | D-2.5 | 705869 | 297.00 |
| ● | 0.5–10 | D-10 | 705870 | 297.00 |
| ● | 2–20 | D-20 | 705872 | 297.00 |
| ● | 5–50 | D-50 | 705873 | 297.00 |
| ● | 10–100 | D-100 | 705874 | 297.00 |
| ● | 20–200 | D-200 | 705878 | 297.00 |
| ● | 100–1000 | D-1000 | 705880 | 297.00 |
| ● | 500–5000 | D-5000 | 705882 | 297.00 |
| ● | 1000–10000 | D-10000 | 705884 | 297.00 |

Transferpette® S-8 and -12 multi-channel pipettes

| | Capacity, μ l (color-coded) | Description | 8-channel Cat. No. | 2021 Price | 12-channel Cat. No. | 2021 Price |
|---|------------------------------------|-------------|-----------------------|---------------|------------------------|---------------|
| ● | 0.5–10 | M8-10 | 705900 | \$760.00 | 705920 | \$940.00 |
| ● | 5–50 | M8-50 | 705906 | 760.00 | 705926 | 940.00 |
| ● | 10–100 | M8-100 | 705908 | 760.00 | 705928 | 940.00 |
| ● | 20–200 | M8-200 | 705910 | 760.00 | 705930 | 940.00 |
| ● | 30–300 | M8-300 | 705912 | 760.00 | 705932 | 940.00 |

Pipette sets and packages, Transferpette® S, adjustable volume, includes tips and shelf/rack (w/sets), benchtop rack (w/packages)



| Type | Description | Cat. No. | 2021 Price |
|--------------------------|------------------------------------|----------|------------|
| Micro | D-2.5, D-10, D-100 | 705890 | \$765.00 |
| Mid | D-20, D-200, D-1000 | 705891 | 765.00 |
| Macro | D-1000, D-5000, D-10000 | 705892 | 765.00 |
| Standard | D-10, D-100, D-1000 | 705893 | 765.00 |
| Pipette package 1 | D-2.5, D-10, D-100, D-200, D-1000 | 705896 | 1,220.00 |
| Pipette package 2 | D-10, D-100, D-200, D-1000, D-5000 | 705897 | 1,220.00 |



Additional product information
can be found at brandtech.com

BRANDTECH[®]
SCIENTIFIC, INC.

888-522-2726
www.brandtech.com