

# VHC<sup>PRO</sup> VacuuHandControl

## VHC<sup>PRO</sup> VacuuHandControl for Fluid Aspiration

Safe, convenient control over your cell culture fluid aspiration. Designed for the VACUUBRAND® BVC series of fluid aspiration systems, the VHC<sup>PRO</sup> also serves as a performance upgrade for other aspiration units.

- Fluid only flows through tip and tubing for leak-proof performance, and easy decontamination.
- Smooth control lever can be held in a number of comfortable positions, in the palm, like a pencil, etc., adapting to your preferred grip.
- Continuous flow setting for prolonged aspiration.
- Works with standard 200 $\mu$ L pipette tips, Pasteur pipette, and serological pipette.



The VHC<sup>PRO</sup> has a number of useful accessories – all autoclavable at 121°C



**8-Channel Manifold with Tip Ejector** features individual removable tip cones for easy maintenance and cleaning, as well as O-rings for low force tip mounting and ejection. Fits standard 200 $\mu$ L tips (not included). Includes multichannel stand, reagent reservoir, wrench and instructions.



**Table Stand** frees up your hands for other work, while keeping tip free from contamination. Smooth surfaces allow for easy wipe down and decontamination.



**Tip Ejector Attachment** allows contact-free removal of 200 $\mu$ L pipette tips to help reduce the chance of cross-contamination (pipette tips not included).

Description	BrandTech Cat. No.	2018 List Price
VHC <sup>PRO</sup> VacuuHandControl Aspiration Controller	688061	\$245.00
Tip Ejector Attachment for VHC <sup>PRO</sup>	635638	79.40
Table Stand for VHC <sup>PRO</sup>	635680	151.00
8-Channel Manifold with Tip Ejector for VHC <sup>PRO</sup>	635679	434.00
Spare silicone tubing for VHC <sup>PRO</sup> (minimum order 2 meter)	636156	20.72/meter

**BRANDTECH**<sup>®</sup>  
SCIENTIFIC, INC.

Toll Free: 888-522-2726  
[www.brandtech.com](http://www.brandtech.com)