

# BRAND® Transferpette® electronic

## Single & Multichannel Pipettes

German manufactured; unprecedented quality and value in an electronic pipette



World's first independently certified "Ergonomics Approved" single and multichannel pipettes!



### Intuitive Interface

User-friendly design allows easy operation, even for new lab techs.

### Adjustable Finger Rest

The adjustable finger rest adapts the instrument to your most comfortable hand position.

### Rechargeable Battery

Over 4000 pipetting cycles between charges means well over a full day of uninterrupted pipetting.

### Low strain

Improved O-rings and the stepped design of the tip ejector reduce tip mounting and ejection forces. O-rings support broad tip compatibility.



### Easy Maintenance

Multichannel manifolds disassemble easily for cleaning or service. See page 6 for more details on the manifolds of BRAND®S mechanical and electronic pipettes.

### Transferpette® electronic single channel pipettes

Capacity, $\mu$ L	without charger		with US plug charger	
	Cat. No.	2018 List Price	Cat. No.	2018 List Price
0.5-10	705339	\$479.00	705319	\$510.00
2-20	705340	479.00	705320	510.00
20-200	705343	479.00	705323	510.00
100-1,000	705346	479.00	705326	510.00
500-5,000 (0.5-5mL)	705347	479.00	705327	510.00

### Transferpette®-8 and -12 electronic multichannel pipettes

Capacity, $\mu$ L	8-channel		12-channel	
	Cat. No.	2018 List Price	Cat. No.	2018 List Price
0.5-10	705419	\$917.00	705469	\$1,097.00
1-20	705420	917.00	705470	1,097.00
5-100	705423	917.00	706573	1,097.00
10-200	705424	917.00	705474	1,097.00
15-300	705426	917.00	705476	1,097.00

**BRANDTECH®**  
SCIENTIFIC, INC.

Toll Free: 888-522-2726

[www.brandtech.com](http://www.brandtech.com)