



Transferpette® S



Transferpette® electronic



Transferpette® S-8/-12



Transferpette® -8/-12 electronic



| Volume range | Transferpette® single channel Nominal volume | | | | | | | | | | | | | Transferpette® multichannel Nominal volume | | | | |
|--------------------|---|--------|-------|-------|----------|-------|-------|--------|--------|--------|---------|------|-------|---|-------|--------|--------|--------|
| | 1 µL | 2.5 µL | 10 µL | 20 µL | 20 µL*** | 25 µL | 50 µL | 100 µL | 200 µL | 500 µL | 1000 µL | 5 mL | 10 mL | 10 µL | 50 µL | 100 µL | 200 µL | 300 µL |
| 0.1 - 20 µL (A) | ✓ | ✓* | ✓* | ✓* | | | | | | | | | | ✓ | | | | |
| 0.5 - 20 µL (B) | | ✓ | ✓ | ✓ | | | | | | | | | | ✓ | | | | |
| 1 - 50 µL (C) | | ✓* | ✓ | ✓ | | | | | | | | | | ✓ | | | | |
| 2 - 200 µL** (D) | | | | | ✓ | ✓ | ✓ | ✓ | ✓ | | | | | | ✓ | ✓ | ✓ | ✓ |
| 5 - 300 µL (E) | | | | | ✓ | ✓ | ✓ | ✓ | ✓ | | | | | | ✓ | ✓ | ✓ | ✓ |
| 50 - 1000 µL** (F) | | | | | | | | | | ✓ | ✓ | | | | | | | |
| 0.5 - 5 mL (G) | | | | | | | | | | | | ✓ | | | | | | |
| 1 - 10 mL (H) | | | | | | | | | | | | | ✓ | | | | | |
| 0.1 - 1 µL (I) | ✓ | ✓* | ✓* | | | | | | | | | | | ✓ | | | | |
| 0.5 - 10 µL (J) | | ✓ | ✓ | ✓ | | | | | | | | | | ✓ | | | | |
| 1 - 20 µL (K) | | | ✓ | ✓ | | | | | | | | | | ✓ | | | | |
| 2 - 20 µL (L) | | | | | ✓ | ✓ | ✓ | ✓ | ✓ | | | | | | ✓ | ✓ | ✓ | ✓ |
| 5 - 100 µL (M) | | | | | ✓ | ✓ | ✓ | ✓ | ✓ | | | | | | ✓ | ✓ | ✓ | ✓ |
| 5 - 200 µL (N) | | | | | | | | | ✓ | | | | | | ✓ | ✓ | ✓ | ✓ |
| 50 - 1000 µL (O) | | | | | | | | | | ✓ | ✓ | | | | | | | |

✓ = Tip volume less than pipettes nominal volume
 * = Calibration and possible adjustment needed
 ** = Bulk tips are yellow or blue colored. racked tips are colorless in a yellow-colored or blue-colored tip tray
 *** = Transferpette® with yellow color code